20421 • TransWarmer® Infant Transport Mattress with WarmGel®

Directions for Use (English)

INTENDED USE: Provides warmth during transport of infant within the hospital or between hospitals.

WARNING: Monitor infant core temperature and check infant's skin for redness. Consult a physician when using this product on premature infants or in non-transport environments.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

DIRECTIONS FOR USE

Do Not Reuse • Disposable • Store at room temperature

Temperature when activated: 40 °C (from 24 °C start), 104 °F (from 75 °F start)

To Activate:

1. Locate metal disc and grasp with thumb and forefingers of both hands.
2. Flex (bend) disc rapidly until crystals begin to form.
3. Massage pack to soften and to increase the activation rate.
4. After pack is fully activated, indent center with palm to form a small nest.
5. Lay infant in nest on the non-woven fabric surface.
6. Check infant skin occasionally for redness and monitor infant core temperature.

CONTENTS

- Sodium acetate (food grade), water and thickeners.
- Do not swallow.
- All contents are non-toxic; however, if contents come in contact with skin or eyes, flush with water.

20421 • TransWarmer® Infant Transport Mattress with WarmGel®

Directions for Use • English (Continued)

EXPLANATION OF SYMBOLS

Reorder Number

Do Not Reuse

Batch Code

ATTENTION: See instructions for use

Latex Free

Product conforms to the Medical Device Directive 93/42/EEC

Authorized Representative in the European Community.

MR Conditional

The **TransWarmer® Infant Transport Mattress** was determined to be MR conditional.

Non clinical testing demonstrated that the TransWarmer Infant Transport Mattress is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

Static magnetic field of 3 Tesla or less

Maximum spatial gradient magnetic field of 720 Gauss/cm or less

MRI-Related Heating

In non clinical testing, the TransWarmer Infant Transport Mattress produced the following temperature rise during MRI performed for 15 min of scanning (i.e., per pulse sequence) in the 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR system:

Highest temperature change +1.6°C (2.88 °F)

Therefore, the MRI related heating experiments for the TransWarmer Infant Transport Mattress at 3 Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7 W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.6°C (2.88 °F).

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the TransWarmer Infant Transport Mattress. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 35 mm relative to the size and shape of the TransWarmer Infant Transport Mattress.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	4,049-mm ²	3,600-mm ²	6,031-mm ²	4,616-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

TransWarmer® and WarmGel® are registered trademarks of CooperSurgical, Inc.

© 2013 CooperSurgical, Inc.

Manufactured for:

CooperSurgical

Made in the USA

95 Corporate Drive
Trumbull, CT 06611 USA
Phone: (800) 243 2974
Fax: (800) 262 0105
www.coopersurgical.com

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TransWarmer® e WarmGel® sono marchi registrati di CooperSurgical, Inc.