


CAUTIONS:

1. Inspect the Spider-Strap™ before each use.
2. Discard the strap if you observe: broken stitching, frayed or cut fabric or straps, cracked or deteriorated hook and loop and/or other signs of deterioration.
3. Do not overtighten to assure adequate circulation.
4. NOT INTENDED FOR VERTICAL LIFTING.
5. DO NOT PLACE CHEST STRAP OVER PATIENT'S DIAPHRAGM.



Federal law restricts this device for sale by or on the order of a physician or other licensed health care professional



Explanation of Symbols	
Symbol	Explanation
	Manufacturer
	European Authorized Representative
	Refer to for Additional Information
	Caution, Consult Accompanying Documents



EMERGO EUROPE
Molensstraat 15
2513 BH, The Hague
The Netherlands
Tel: +31 (0)70 345 8570
Fax: +31 (0)70 346 7299



Emergency Products & Research
890 West Main Street
Kent, OH 44240-2218
+1-330-673-5003

CAT #8-05
or www.epandr.com

Manufactured Exclusively By



Creating Industry Standards.

Emergency Products & Research, Inc.
890 West Main Street
Kent, OH 44240-2218

1-800-322-5725
www.epandr.com

©2006 Emergency Products & Research, Inc.

Life-Assist
Helping Heroes Save Lives

SPIDER-STRAP™
Immobilization Strap



User Manual (EP-499)

(Spine Board shown for illustration purposes only, not included with Spider Strap™)


Manufactured Exclusively By



Creating Industry Standards.

SPIDER-STRAP™

Immobilization Strap



V-Strap
(over and
down on the
shoulders)

Chest Strap
(Do not place
over diaphragm,
up to axilla)

Pelvic Strap
(Across hip bones)

Thigh Strap
(High on thighs)

Lower Leg Strap
(Below the knee)

For best results in the field,
practice using the Spider-Strap™
until you become thoroughly
familiar with it.

- ### APPLICATION:
1. Hold folded strap above patient's chest with strap loops at bottom pointing to patient's feet.
 2. Unfasten small retaining fastener which holds entire strapping system together.
 3. Using both hands, divide the strap in the middle, letting it fall open naturally.
 4. Place strap on patient's chest with the center strap folds pointing towards the patient's feet and the V-Strap on the patient's chest.
 5. Pull the top cross strap (foot strap) towards the patient's feet. If it has been done correctly, you will see a solid color strap down the center of the patient. If you see any other colored straps crossing the center strap, that strap is inverted.
 6. Unfasten the hook and loop on each cross strap and run through the holes in the backboard or scoop stretcher and then back around and fasten it to itself.

- ### REFOLDING:
1. Lay strap out as if on patient.
 2. Fold cross straps back on themselves. Secure with hook and loop.
 3. Position yourself at the foot end of the strap and grasp the right and left sides of the foot strap.
 4. Pick up cross straps in sequence, half in your right hand and half in your left hand, beginning with the bottom strap. Align the center strap folds toward you.
 5. When entire strapping system has been picked up, fold in half and secure with the small hook and loop fastener located on the V-Strap.
 6. Store near backboard or scoop stretcher until next use.
- ### CARE:
1. Keep hook and loop clean; the hook and loop closures will support considerable weight and "safe use" if kept clean and lint free.
 2. Hand wash. For machine washing, put strap in a pillow case, use cool water and a mild detergent and the gentle cycle.
- For Professional Use Only!**

**AIR DRY ONLY. DO NOT PUT
IN DRYER. DO NOT AUTOCLAVE.**

Manufactured Exclusively By
EMERGENCY
Products + Research
Creating Industry Standards.

Life-Assist
Helping Heroes Save Lives