

Ethox

INFU-SURG® PRESSURE INFUSION & INFU-SURG® CLARITY BAG 4005H, 4010H, 4030, 4050 & 4100

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK



FIGURE 1: Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use



FIGURE 2: Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



FIGURE 3: 360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation

INSTRUCTIONS FOR USE

1. Remove pressure infusion bag from package, unfold pressure bag and hang pressure bag on to IV pole by the loop located at the top of the pressure bag (Figure 1).
2. Insert intravenous bag between mesh panel and pressure bag on 4050H and 4010H only.
3. If using clarity bag 4005 or 4100, place intravenous bag between clear panel and pressure bag.
4. Hang loop of intravenous bag on to welded solution bag hook located just above panel on pressure bag (Figure 1).
5. FOR 4030 ONLY: Slide hanger loop through top of 3000 ml intravenous bag and rehang on IV pole.
6. Inflate pressure infusion bag by rotating the off portion of the valve/stopcock to the right. Pump inflation bulb until gauge reaches the proper pressure (300 mmHg) indicated in the green portion on the manometer (Figure 2).
7. Once proper pressure is reached on gauge, rotate off portion of the valve/stopcock straight up towards the pressure bag to maintain pressure in bag. **NOTE: precise positioning of valve/stopcock is important or bag could lose pressure.**
8. To deflate bag turn the off portion on the valve/stopcock straight down towards inflation bulb to release pressure in bag.
9. Proper valve positioning is illustrated on all pressure bags.

All Infu-Surg Pressure Infusion Bags are recommended for single patient use.

CAUTION:

All Infu-Surg Pressure Infusion bags come equipped with a 330 mmHg pressure relief valve that prevents over inflation. Once 330 mmHg is reached, the manometer will indicate red on the manometer along with the word "caution" and pressure will begin to release from valve. Pressure will continue to release from valve until pressure decreases below 330 mmHg (Figure 3).