# Stifneck Select

Select and Pedi-Select Extrication Collar

Directions for Use

- 980010 Stifneck Select (Qty.1)  
980012 Stifneck Select, green (Qty.1)  
980020 Stifneck Pedi-Select (Qty.1)  
989910 Stifneck Select, Private Label (Qty.1)  
989920 Stifneck Pedi-Select, Private Label (Qty.1)

**Manufactured in China for:**  
LAERDAL MEDICAL AS, Norway  
P.O. Box 377, N-4002 Stavanger, Norway  
Tel. +47 51 51 17 00,  
Fax +47 51 52 35 57

**Distributed by:**  
Canada and Latin America:  
LAERDAL MEDICAL CANADA LTD.  
151 Nashdene Rd, Unit #45  
Toronto, ON, Canada, M1V 4C3  
Tel. +1 (416) 298-9600,  
Toll free 888/LAERDAL (523-7325)  
ou en français (800) 567-9987  
Fax +1 (416) 298-8016

**USA:**  
LAERDAL MEDICAL CORPORATION  
167 Myers Corners Road, P.O. Box 1840  
Wappingers Falls, New York 12590-8840  
Tel. (800) 431-1055, +1 (845) 297-7770  
Fax (800) 227-1143, +1 (845) 298-4545

**Australia:**  
Laerdal Pty Ltd.  
8 Stamford Road  
Oakleigh VIC 3166.

**Laerdal do Brasil**  
Al.Tocantins 125 – sl 01 – Alphaville  
Barueri / SP / Brazil  
06543-260  
Tel: +55 11 4193-8007  
e-mail: comercial@laerdal.com

[www.laerdal.com](http://www.laerdal.com)

ENGLISH  
ČESKÝ  
DANSK  
DEUTSCH  
ESPAÑOL  
FRANÇAIS  
ITALIANO  
NEDERLANDS  
NORSK  
POLSKI  
PORTUGUÊS  
Русский  
SVENSK  
SUOMI  
日本語  
中文  
한국어

**Laerdal**  
helping save lives

## Contents

---


ENGLISH	4
ČESKÝ	5
DANSK	6
DEUTSCH	7
ESPAÑOL	8
FRANÇAIS	9
ITALIANO	10
NEDERLANDS	11
NORSK	12
POLSKI	13
PORTUGUÊS	14
Русский	15
SVENSK	16
SUOMI	17
日本語	18
中文	19
한국어	20
REGULATORY INFORMATION	21

# English

## Intended Use

This product is intended to support cervical spine in a neutral position during transportation, in combination with other cervical and full body immobilization devices (CSIDs and FBSIDs).

Cannula hooks are provided for the securement of supplemental oxygen lines.

## Important


- In cases of a suspected spinal cord injury, proper neck immobilization is only one part of total immobilization.
- It is imperative that the patient be properly immobilized to prevent any movement of the spinal column.

## Warnings

- Stifneck Select Collars should only be used by qualified healthcare providers.
- Improper use may result in injury or permanent disability. Always follow local spinal immobilization protocol.
- Use the correct size collar. An oversized collar may hyperextend a patient's cervical spine. An undersized collar may not provide appropriate stability.
- Do not adjust the Select or Pedi-Select collar on patient.

## User Instructions


### 1. Measure the Patient


Align the head to neutral or "eyes forward" position unless contraindicated by local protocol.

### 2. Match the collar size to the patient:


Select: Choose from 4 adult positions.  
Pedi-Select: Choose from 3 Paediatric positions.


### 3. Adjust and Lock the adjustable collar


Use the same method for adjusting both Collars.

### 4. Preform the collar


### 5. Apply the collar while manually maintaining neutral head position

Place the chin support under the chin. If a different size is needed, remove, resize, and re-apply the collar.

Pull the back of the collar snug while holding the front in place, and then fasten.

For a supine patient, slide the rear panel behind the neck before placing the chin support.

## Technical Specifications

### Storage:

Do not store collar in folded position.

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© 2015 Laerdal Medical AS. All rights reserved