ENGLISH

Brand Name: Laerdal BaXstrap Spineboard
Common Name: Backboard/Spineboard/Longboard
Description: Full-body length spinal immobilization device

Intended Use: A base structure to be used with other adjunct cervical spine and head immobilization devices to facilitate in-line, neutral immobilization and transport of adult and pediatric patients.

Indications for Use: For use with patients having known or suspected spine or spinal cord injuries when patient immobilization and transportation is indicated.

Contraindications: None known

- Notes
- Always follow local protocol for patient movement and immobilization techniques.
 - Not made with natural rubber latex.
 - The product is X-ray, CT, and MRI compatible.

Warnings and Cautions: Read and thoroughly understand these directions. These directions are not intended as, and do not in any way constitute medical treatment advice, instruction in or teaching of medical practice, or imply or convey any rights, medical direction or authority. This and other spinal immobilization adjunct devices should only be used by personnel who have received adequate training. Inspect the BaXstrap prior to initial use and before each patient use. Pay particular attention to the underside of the head end. If the product develops punctures, broken pins or cracks then discontinue use immediately.

Directions for Use:

- Immobilize the cervical spine using a proper Cervical Immobilization Device

- Position the patient on the BaXstrap using appropriate number of personnel.

- Properly fasten/strap the patient to the BaXstrap to protect against patient movement in any direction (up, down, left or right).

- Secure the head to the BaXstrap using a proper Head Immobilization Device.

- Assure in-line neutral immobilization of the head, neck, torso and pelvis.

Cleaning: Non-sterile, non-critical infection risk. Withstands cleaning in soap and water (1:16) and 1:100 household bleach: water.

Environmental Conditions:

Temperature (Operating and Storage):

-34 °C to +52 °C (-29 °F to +126 °F)

Humidity (Operating and Storage):

5-95% RH non-condensing.

The product has a 5 year limited warranty. See the enclosed "Laerdal Global Warranty" for terms and conditions. The warranty is also available at www.laerdal.com.

 This medical device complies with the general safety and performance requirements of Regulation (EU) 2017/745 for medical devices.

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 Medicinteknisk produkt

© 2020 Laerdal Medical AS.All rights reserved.

US Patent 5,950,627
BaXstrap is a trademark of Laerdal Medical AS.
EN 1865 compliant

Manufacturer: in USA for Laerdal Medical AS

P.O. Box 377
Tanke Svländsgate 30, 4002 Stavanger, Norway
T: (+47) 51 51 17 00

Printed in USA

Distributed in USA by:
LAERDAL MEDICAL CORPORATION
167 Myers Corners Road,
Wappingers Falls New York 12590
Tel. +1 (800) 431-1055, +1 (845) 297-7770
Fax +1 (800) 227-1143
E-mail: customerservice@laerdal.com

Distributed in Canada by:
LAERDAL MEDICAL CANADA LTD.
151 Nashdene Rd., Unit #45,
Toronto, ON, Canada, M1V 4C3
Tel. +1 (416) 298-9600
Fax +1 (416) 298-8016
E-mail: savelives@laerdal.ca

生产厂家:挪威 Laerdal Medical AS

注册地址:
P.O. Box 377, Tanke Svländsgate 30,
4002 Stavanger, Norway

生产地址:
Formed Plastics Inc, 207 Stonehinge Lane Carle
Place, NY 11514, USA

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5585 Rev. F

www.laerdal.com  **Laerdal**
helping save lives

Life-Assist
Helping Heroes Save Lives

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