FAST. ACCURATE. EASY TO USE.

microdot[®] Xtra

Blood Glucose Monitoring System

OPERATIONS & QUALITY ASSURANCE PROCEDURE MANUAL




TABLE OF CONTENTS

Chapter	Page
1. microdot® Xtra Meter Overview & System Components	2 - 4
2. Initial Setup of the microdot® Xtra Meter	5 - 7
3. Control Solution Testing	8 - 10
4. Testing with Patients	11 - 12
5. microdot® Xtra Meter Specifications	13
6. In-Service Training Outcomes	14
7. microdot® Xtra Blood Glucose Meter Log Sheets.	
7A. microdot® Xtra BGMS In-Service Form	_
7B. microdot® Xtra BGMS Qualified Trainers Form	17
7C. microdot® Xtra BGMS Qualified Trainer/Operator	
Certification Form	
7D. microdot® Xtra BGMS Quality Control Record	19
8. Material Safety Data Sheets	
8A. microdot® Xtra Test Strips	20
8B. microdot® Control Solutions	
8C. microdot® Bleach Wipe	25 - 28
8D. microdot® Minute Wipe	29 -36
9. Troubleshooting the microdot® Xtra Meter	37
10. Manufacturer Cleaning Instructions	38 - 40

microdot° Xtra BLOOD GLUCOSE METER OVERVIEW

This section provides general information on the microdot® Xtra Blood Glucose Monitoring System

The microdot® Xtra Blood Glucose Monitoring System is a portable whole blood testing system which performs tests with quick and accurate results.

It consists of three main parts:

- 1. microdot® Xtra Blood Glucose Meter
- 2. microdot® Xtra Test Strips
- 3. microdot® High and Low Control Solution

These products have been designed, tested and proven to work together as a system to produce safe and accurate results.

<u>Use only microdot® Xtra Test Strips</u> <u>and Control Solutions with the</u> <u>microdot® Xtra Glucose Meter.</u>

Your system includes:

- microdot® Xtra Meter
- Protective Orange Meter Sleeve
- Quality Assurance Manual
- Lifetime Warranty
- 3V Lithium Battery (installed)

Additional components include:

- microdot® Xtra Test Strips
- microdot® High and Low Control Solutions

SYSTEM COMPONENTS


microdot®Xtra Blood Glucose Meter

LCD Screen: Shows blood glucose result and symbols that guide you through the test.


Up/Down Toggle Buttons: Scrolls the memory, sets time and date.

Eject Button: Releases test strip.

Test Strip Port: Insertion site for test strip. -



Battery Compartment: Holds one 3V Lithium battery (CR2032).



Customer Service Number: Located on the back of the meter.

Call this number if you have any questions or problems with the microdot® Xtra Glucose Meter.

(877) 374-4062

SYSTEM COMPONENTS

microdot® Xtra Test Strip


The microdot® Xtra Blood Glucose Test Strip offers the latest advances in biosensor, auto-code technology. Blood is applied to the top edge of the microdot® Xtra Test Strip and is automatically drawn into the white channel where the reaction takes place. Only 600 nanoliters (0.6 microliters) is needed.

The microdot® Xtra Test Strip consists of the following parts:

Top Edge: Apply a drop of blood here, where the white channel meets the top of the edge strip.

White Channel: This is where you check if enough blood has been applied to the top edge.

Contact Bars: Insert this end of the test strip into the meter. Push firmly until the strip can go no further.




INITIAL SETUP OF THE microdot® Xtra METER

Before using the microdot® Xtra Meter for the first time, you should set the actual Time, Date and Year. The unit of measurement is preset and cannot be changed.

Setting the Time, Date and Year

Enter the Set Mode.

To Enter the set mode, turn the meter on by pressing the **C button**. After the segment test, the time and date will start to flash.



Step 1

Set the Time and Date Format.


Pressing the **C button** again will now toggle between U.S. and International Time and Date formats.



For U.S.: 12h Time format, mm-dd (begins with AM setting). To accept the desired setting, press the **M button**.


Step 2 Set the Hour.

The Hour will start to flash. It can now be changed by pressing the **C button**. The U.S. setting will begin at 12:01 AM and will move to PM. To accept the correct setting, press the **M button**.



Step 3 Set the Minutes.

The Minute will start to flash, it can now be changed by pressing the **C button**. To accept the correct setting, press the **M button**.



Chapter 2

INITIAL SETUP OF THE microdot® Xtra METER

Setting the Beeper

After setting the Time and Date, the Sound Symbol will appear and can now be changed by pressing the C button. To accept the setting, press the M button.



ON

This option is used to switch the beeper on or off . When turned ON, a sound will be heard when blood or control solution is applied to the strip and when the test is finished. The sound will be heard when an error has occurred or if an alarm is triggered.



When the M button is pressed to accept the beeper option, the meter will display END and switch off.



INITIAL SETUP OF THE microdot® Xtra METER

Using Meter Memory


Your microdot® Xtra Meter stores the 500 most recent blood glucose, control solution test results and insulin data with date and time in the memory. It also provides you with 14-day averages of your blood glucose test results. You can review the test results in memory with these easy steps.

Step 1 Enter the Memory Mode.

To Enter Memory Mode, turn on the meter by pressing the **M button**.

The meter will display the last result with Mem. Symbol, Time and Date.

1.






Step 2 Recall Test Results.

Previous results can be displayed by pressing the **C button**. As long as the **C button** is held, the meter will scroll through the memory displaying the result and the memory location. Once the button is released, the meter will display the selected result with its Time and Date. When the memory is full (500 results stored), the oldest result is dropped as the newest is added.

NOTE: When using the meter for the first time "Mem.---" will appear, showing that there are no test results stored in memory.

2.





CONTROL SOLUTION TESTING

This chapter describes the necessary steps to test with the microdot[®] Control Solutions in order to validate the performance of the microdot[®] Xtra Meter and microdot[®] Xtra Test Strips.

When Should You Conduct a Control Solution Test?

- Any time you open a new vial of test strips.
- Whenever you think the system is not working properly.
- If the blood glucose test results differ from the resident's symptoms or non-symptoms.
- If you believe the results are not accurate.
- If you drop the meter.
- If the vial of test strips has been left open for an extended period of time.

It is critical to follow the Operating Guidelines, on the right, to obtain accurate results while using the microdot® Xtra Glucose System.

Operating Guidelines:

- Use only microdot® Control Solutions (High and Low).
- Check the expiration date on the control solution vial. Do not use if expired or if the discard date has passed.
- Control solution, meter and test strips should come to room temperature before testing (66-77°F / 20-25°C).
- Use solution for three months after first opening. Record the discard date (opening date plus three months) on the control solution vial. Discard after three months.
- Close tightly and store the control solution at temperatures between 50-86°F (10-30°C).
- Do not refrigerate. Do not freeze.


CONTROL SOLUTION TESTING

Notes:


The control solution test is similar to a blood test except that you use microdot® Control Solution instead of a drop of blood. The control solution ranges printed on the test strip vial are for microdot® Control Solutions only. It is used to check the meter and test strip performance. It is NOT a recommended range for your blood glucose level. When performing a control solution test it does not matter which solution you use first, the Low or High.

Control Solution Test Procedure:

- **1.** Shake Low Control Solution bottle well before using.
- 2. Remove cap and discard the first drop of Control Solution and wipe off the dispenser tip to ensure a good sample and an accurate result.
- **3.** Insert a test strip into the microdot® Xtra Meter. Be sure the black contact bars go into the meter. Push the strip in firmly. Be sure it can go no further.



4. Invert bottle and squeeze out one drop of control solution. Apply the drop to the strip by bringing the meter and the strip to the drop. Touch the drop with the top edge of the test strip and wait until the test pad fills with the solution. Results appear in 10 seconds.



CONTROL SOLUTION TESTING

- **5.** Compare the results with the ranges of expected results shown on the test strip vial. (Low = Blue Cap, High = Red Cap)
- **6.** You should obtain results within the expected range printed on the test strip vial. If this is not so, repeat the test.

If the results are out of the range printed on the test strip vial, check the following:

- Was the vial at room temperature?
- Did you shake the bottle of control solution before using?
- Has the control solution expiration or discard date expired?
- *Is the meter malfunctioning?*
- **7.** Repeat steps 1-6 for the High Control Solution procedure.

If the test result is still out of range, call your microdot® Customer Service Representative at:

Toll Free (877) 374-4062

DO NOT test blood until you obtain control results within the expected ranges.

TESTING WITH PATIENTS

This section describes the procedure to test patient blood samples using the microdot® Xtra Blood Glucose System.

Operating Guidelines

- Before attempting to test with real blood, make sure you have performed control solution tests correctly to insure the meter and test strips are performing properly and to verify technique.
- Use the test strips before their expiration date and within three months after opening.
- Do not use test strips that are wet, bent, scratched or damaged. Use each test strip immediately after removing it from the vial.

Patient Test Procedure:

	What You Do	Notes
Mem. Avg. CODE mgdil mmalit mmalit	1. Press the C button to turn the meter on and/or insert test strip by pushing the two contact bars in firmly.	To verify all display symbols are working, all symbols should appear at the same time. If any of the symbols are missing or do not display completely, contact Cambridge Sensors USA at (877) 374-4062
	2. Obtain a drop of blood.	Using a safety lancet, lance the side of the finger to obtain a rounded blood sample. Avoid squeezing the puncture site excessively.

TESTING WITH PATIENTS

	What You Do	Notes					
OK	3. Apply sample.	The blood drop symbol will flash alternately to indicate the meter has detected that a test strip has been inserte and is ready for the blood sample to be added.					
		Apply the drop of blood directly to the edge of the test strip, and allow the sample to wick automatically into the test zone.					
NOT OK		If you do not apply a blood sample within one minute, the meter will turn itself off. Either reinsert the test strip or press the M button to turn the meter back on.					
108 mg/dL	4. Results in 10 seconds.	After an adequate blood sample has been applied to the test strip, the microdot® Xtra Meter will display three running dashes, which indicates the meter is performing the test. After 10 seconds, the test result will display.					
	5. Dispose of lancet.	Dispose of used lancet in an approved sharps container.					
	6. Dispose of test strip.	Simply press the eject button on the meter to dispose of used test strip.					
Note: If the blood test result "Lo" will appear on the mete This indicates severe hypogl (low blood glucose). You sho immediately treat your hypoas recommended by your ph	er display. ycemia ould oglycemia	If the blood test result is higher than 525 mg/dL, "Hi" will appear on the meter display. This indicates severe hyperglycemia (high blood glucose). You should immediately treat your hyperglycemia as recommended by your physician's					

protocol.

protocol.

microdot® Xtra METER SPECIFICATIONS

0.6 Microliters (600 Nanoliters) Sample Size

10 Seconds **Total Test Time**

Strip is placed in the meter and sample **Sample Application**

wicks into strip

Hematocrit Range 30 - 50%

Measurement Range 20 to 525 mg/dL

Interference No interference from over 20 interfering substances

Strip Removal Strip ejection by push button

Coding Auto-code

Accuracy of Strip/Meter System At least +/- 20% relative to YSI in clinical trials

Correction Coefficient of Regression > 0.969

Precision of Strip The strip variation is not greater than 6.4%

Shelf Life of Meter Approx. 5 years

10 - 90%; storage 90% max (non-condensing) **Operating Humidity range**

50° - 104°F (10° - 40°C) **Operating Temperature** Plasma equivalent **Type of Glucose Result**

Strip Vial Packing Size 50 strips per vial

Typical Control Solution Ranges

Printed on vial label when used with Strips/Meter

3 months **Strip Shelf Life after Opening**

Strip Vial Plastic with desiccant sleeve vial

Altitude To 10,000 ft (Target)

Module Size 65 mm x 85 mm x 16 mm

Approx. 60 grams / 2 ounces Weight

Battery Life 1,500 tests or about 1 year at three tests per day

Glucose Units mg/dL or mmol/L

Total of 500 glucose results and insulin-input data Memory

Power Source One replaceable 3V lithium battery

(CR2032 or equivalent)

One minute after last user action **Automatic Shut-off**

1 Year Manufacturer Warranty Warranty

IN-SERVICE TRAINING OUTCOMES

Once your microdot® Xtra Meter in-service is complete, Health Care Professionals should be able to:

1. Locate and explain the following components of the microdot® Xtra Meter:

> Eject Button Test Strip Port **Battery Compartment** Serial Number **Customer Service Number**

2. Locate and explain the label information on the microdot® Xtra Test Strip Vial:

> Discard Date Lot Number **Expiration Date Control Ranges** Storage Temperature Range

- **3.** Set Time and Date of meter
- **4.** Identify three parts of the microdot® Xtra Test Strips

Top Edge Test Pad **Contact Bars** **5.** Identify the following parts of the microdot® Control Solution Bottles:

> Low vs. High **Discard Date** Lot Number **Expiration Date** Storage Temperature Range

- 6. Properly insert the microdot® Xtra Test Strip
- 7. Perform control solution tests
- 8. Document and maintain microdot® Xtra **Quality Control Records**
- 9. Obtain a blood sample
- **10.** Perform a blood test
- **11.** Identify meter result range
- **12.** Access meter memory
- **13.** Change the battery
- **14.** Identify and resolve error readings

microdot® Xtra BLOOD GLUCOSE MONITORING SYSTEM **LOG SHEETS**

The forms in this chapter are Master Copies. Please make photocopies for distribution. If you need replacement documents, please contact Customer Service at (877) 374-4062

microdot®xtra	BLOOD GLUCOSE MONITORING	SYSTEM
IN-SERVICE FORM		

The Health Care professionals listed below are trained and have demonstrated
proficiency using the microdot® Xtra Blood Glucose Monitoring System.

Cambridge Sensors Re	presentative/Qualified Trainer

Date	Health Care Professional's Name

The chart is designed for use during the initial training on the microdot® Xtra Blood Glucose System by a Cambridge Sensors Representative, In-Service Video or Qualified Trainer.

microdot® Xtra BLOOD GLUCOSE MONITORING SYSTEM QUALIFIED TRAINERS FORM

Date	Health Care Professional's Name	Title
hridae Sensors	Representative	
iwi iage Jelijolj	nepresentative	
huidas Canas	Donyos antativo Signatura	
briage Sensors	Representative Signature	

Date

microdot®Xtra BLOOD GLUCOSE MONITORING SYSTEM **QUALIFIED TRAINER/OPERATOR CERTIFICATION FORM**

I. microdot®	Xtra	Meter -	Locate	the	foll	lowing	g:
--------------	------	---------	--------	-----	------	--------	----

- a. Battery
- b. Serial Number
- c. Toll Free Customer Service Number
- d. Eject Button

II. Identify and explain the following:

- a. Blood Test Procedure
- b. Control Solution Procedure
- c. Troubleshooting

III. Explain proper procedure for:

- a. Dating of microdot® Xtra Test Strips and Control Solution
- b. When to perform a Control Test
- c. Handling Control Solution results that are out of range

IV. Perform and explain the following procedures:

- a. Control Solution testing
- b. Blood testing
- c. Recalling test results
- d. Changing the battery

If all tasks were not properly completed, have the operator repeat the procedure correctly. When all tasks have been properly completed, sign bottom of this form and fill out the proper Qualified Trainer or In-service form.

Cambridge Sensors USA Sales Rep	_ Date	
Qualified Trainer	Date	
Operator	Date	_

$microdot^{\circ}$ Xtra BLOOD GLUCOSE MONITORING SYSTEM **OUALITY CONTROL RECORD**

	Q 0 / 1 1													
	• •	• •	•	• •	•	•	• •	•	•	•	•	•	•	• •
microdot® Xtra Meter Serial Number														
Month/Year														

Date	Station/ Initials		Low Control Range (mg./dl.)	Low Control Result (mg./dl.)	High Control Range (mg./dl.)	High Control Result (mg./dl.)	Natar
Date	initiais	Lot #	(mg./dl.)	(mg./dl.)	(mg./dl.)	(mg./dl.)	Notes:
	<u> </u>				l	l	

microdot®xtra Blood Glucose Test Strips

Section I - Identification

Trade Name: microdot® Xtra Blood Glucose Test Strips

Common Name:Blood glucose test strips for use with microdot®

Xtra Blood Glucose Test System

Supplier: Cambridge Sensors USA

10051 Bode Rd. Plainfield, IL 60585

Emergency Phone Number: (877) 374-4062

Section II - Composition

This product does not present a physical or health hazard under reasonable use or under emergency situations involving a release of only this product. This material is therefore not considered to be a "Hazardous Chemical" as defined by the Federal Occupational Safety and Health Administration in the Hazard Communication Standard (29 CFR 1910.1200) or the equivalent standards generated by state agencies.

ACCORDINGLY NO MATERIAL SAFETY DATA SHEET IS REQUIRED FOR THIS PRODUCT.

MSDSs that represent non-hazardous chemicals are not covered by the HCS. Paragraph 29 CFR 1910.1200 (g) (8) of the standard requires that "the employer shall maintain in the work-place copies of the required MSDSs for each hazardous chemical, and shall insure that they are readily accessible during each work shift to employees when that are in their work area(s)". OSHA does not require or encourage employers to maintain MSDSs for non-hazardous chemicals. Consequently, an employer is free to discard MSDSs for non-hazardous chemicals.

microdot® Control Solutions

Section I - Identification

Trade Name: microdot® Control Solutions

Common Name: Control Solutions for use with microdot® Xtra

Blood Glucose Test System

Supplier: Cambridge Sensors USA

> 10051 Bode Rd. Plainfield, IL 60585

Emergency Phone Number: (877) 374-4062

Section II - Composition / Information on Ingredients

Salt solution **Chemical Family:**

Chemical Name: N/A

Section III - Hazards Identification

May cause irritation if in contact with eye. **Eye Contact:**

Skin contact may cause irritation. **Skin Contact:**

Though not a likely route of occupational exposure, **Ingestion:**

ingestion of the product may cause choking, vomiting

or nausea.

N/A Inhalation: N/A **Chronic Exposure**

Section IV - First Aid Measures

Although no adverse health reactions are expected from **Eye Contact:**

> the normal use of this product, it is recommended to flush eyes with water and seek medical advice whenever there

is a potential injury to the eye.

Skin Contact: If contact with this product leads to reddening, inflammation

or irritation, flush the exposed area with running water. If

irritation persists, get medical attention.

As a precaution, get medical attention if there has been **Ingestion:**

ingestion of this product.

If breathing becomes difficult, remove victim to fresh air Inhalation:

and get medical attention.

Chronic Exposure: N/A

microdot® Control Solutions

Section V - Fire Fighting Measures

Auto-flammability Flash Point (test method) **Extinguishing Media Special Fire Fighting Procedures**

Not determined Use fire extinguishing media appropriate for site collections. Structural fire fighting gear and self-contained breathing apparatus will provide adequate protection if this product is in a fire area.

Fire and Explosion Hazards Not determined **Hazardous Combustion**

Products Upper Explosion Limit (%) Lower Explosion Limit (%)

Thermal decomposition may emit carbon monoxide and

carbon dioxide. Not determined Not determined

Not determined

Section VI - Accidental Release Measures

Spill and Leak Procedures

Use an absorbent material to contain/pick up the spilled solution. Place all contaminated disposals into a suitable container, seal, label and hold for disposal.

Section VII - Handling and Storage

Storage Temperature Handling/Storage **Ventilation Requirements**

Sensitivity to Static Electricity Sensitivity to Mechanical Impact Store vials as directed in the package insert.

Handle and store vials as directed by the package insert.

No special requirements.

Not known Not known

Section VIII - Exposure Controls / Personal Protection

Respiratory Protection

Other Engineering Controls

Ventilation

Protective Gloves

Work Practices

Not required under normal use of this product.

Not required under normal use of this product. Wear appropriate gloves to prevent skin contact. Replace

torn or punctured gloves promptly.

Wear appropriate eye protection to prevent eye contact. **Other Protective Equipment**

Wear appropriate body protection to prevent skin contact.

Eye wash stations and deluge showers.

Good laboratory technique should be used when handling

this product. Observe appropriate chemical hygiene.

Do not place in mouth.

Hygienic Practices Do not eat, drink or smoke while working with this product.

Upon completion of work activities involving this product,

wash hands thoroughly with soap and water.

microdot® Control Solutions

Section IX - Physical and Chemical Properties

Pure Substance or preparation Preparation Physical Form Liquid Appearance/Odor Blue, odorless pH AS is neutral

Odor Threshold Not determined **Melting/Freezing** Not determined **Partition** Not determined **Evaporation Rate** Not determined Vapor Pressure (mmHg) Not determined Vapor Density (air =1) Not determined **Viscosity** Not determined **Volatiles** Not determined **Volatile Organic Compounds** Not determined **Auto-flammability** Not determined **Flash Point** Not determined **Oxidizing Properties** Not determined

Stability

Materials to Avoid Strong bases, strong acids and water reactive materials Thermal decomposition may emit carbon monoxide and **Hazardous Decomposition**

Stable

Products

carbon dioxide.

Section X - Toxicological Information

Route of Entry Ingestion, skin and/or eye contact

Effects of Chronic Exposure Not known **Effects of Acute Exposure** Not known **Special Health Effects** Not known **Target Organs** Not known

Section XI - Ecological Information

Potential Effect on Environment Not known Not known **Potential to Bioaccumulate** Not known **Mobility** Not known **Ecotoxicity Persistence and Degradeability** Not known **Agua Toxicity** Not known

Section XII - Disposal Considerations

Water Disposal Please consult local, state and federal regulations for

additional guidance on disposal.

Empty Container Warnings Not known

Chapter 8

SAFETY DATA SHEET

microdot® Control Solutions

Section XII - Transport Information (See also Section IX)

ADD / DID	
ADR / RID	Not known
CEFIC Tremcard	Not known
Hazchem Code	Not known
Kemmler Code	Not known
IMDG Classification	Not known
IATA Classification	Not known
Marine Pollutant	Not known
UN Number	Not known
UN Class	Not known
UN Packing Group	Not known

Section XIV - Regulatory Information

EEC Hazard Classification	Not known
Risk Phrases	Not known
Safety Phrases	Not known

Section XV - Other Information

Directives 88/379/EEC and 91/155/EEC have been considered when compiling this MSDS; the information is provided for health and safety assessment by an industrial user. Reference should be made to any relevant local or national health, safety or environmental legislation. This information does not constitute indication of suitability for specific uses. The information, data, and recommendations contained herein are based upon information believed by Cambridge Sensors USA, after reasonable investigation and research, to be accurate. However, Cambridge Sensors USA does not warrant the accuracy of this information. All materials and mixtures may present unknown hazards and should be used with caution. When necessary or appropriate, independent opinions regarding the risk of handling or exposure should be obtained from trained professionals. Cambridge Sensors USA disclaims any warranty against patent infringement and the implied warranties of merchantability and fitness for a particular purpose. Customer's sole and exclusive remedy shall be replacement of the product or return of the product and refund of the purchase price, at Cambridge Sensors USA option. In no case can Cambridge Sensors USA be liable for incidental or consequential damages, including lost profits.

microdot® Bleach Wipe

Product and Company Identification

microdot® BLEACH WIPE **PRODUCT NAME:**

PRODUCT DESCRIPTION: WIPE PRESATURATED WITH PRE-DILUTED BLEACH

microdot® BLEACH WIPE TRADE NAME: **GENERAL USE:** Hospital Disinfectant Sodium hypochlorite CHEMICAL FAMILY:

Bleach disinfectant formula absorbed on towels. Mild bleach odor. PRODUCT DESCRIPTION:

microdot® BLEACH WIPE

HOSPITAL CLEANER DISINFECTANT TOWELS WITH BLEACH

MANUFACTURER: MEDTROL, INC.

CAMBRIDGE SENSORS USA, LLC **MANUFACTURED FOR:**

DATE PREPARED: September 25, 2012

ADDRESS: 7157 N. AUSTIN AV. NILES, IL 60714 10051 BODE RD. PLAINFIELD, IL 60585 **ADDRESS:**

TELEPHONE NUMBER FOR

INFORMATION / CUSTOMER SERVICE: 877-374-4062

24-HOUR EMERGENCY

TELEPHONE NUMBER: CHEMTREC 800-424-9300

Section II - Composition/Information on Ingredients

Hazardous Components % (by Weight) CAS# <1.05% 7681-52-9 Sodium hypochlorite

Notes: * The balance of ingredients not listed above are non-hazardous as defined in the OSHA hazard communication standard 29CFR 1910.1200.

Section III - Hazards Identification

EMERGENCY OVERVIEW:

Avoid contact with eyes, skin, and clothing as this product may produce irritation. Do not allow this product to contact acidic materials as hazardous chlorine gas may be released.

POTENTIAL HEALTH EFFECTS

INHALATION: No adverse effects are anticipated from inhalation.

SKIN: Normal exposure (contact) is not likely to cause significant skin irritation.

EYES: Causes moderate eve irritation.

INGESTION: May cause gastrointestinal irritation and upset.

CARCINOGENICITY:

NTP? NO IARC MONOGRAPHS? NO OSHA REGULATED? NO

CALIFORNIA, Prop.65? NO

microdot® Bleach Wipe

Section IV - First Aid Measures

INHALATION: No specific treatment - suspected hazard by this route is minimal. Note to Physician:

No specific antidote. Supportive care. Treatment based on judgment of the physician

in response to reactions of the patient.

EYES: If in eyes, hold eye open and rinse slowly and gently with water for 15-20 minutes.

Remove contact lenses, if present, after the first five minutes, then continue rinsing eyes.

Seek medical attention if irritation persists.

SKIN: Wash with soap and water. No further first aid should be required. **INGESTION:** No specific treatment - suspected hazard by this route is minimal.

Section V- Fire Fighting Measures

GENERAL HAZARDS: Product is minimally flammable FP≥ 101°C (215°F). Chlorine may be released

from this product in the presence of acids.

EXTINGUISHING MEDIA: Water fog, carbon dioxide and dry chemical to fight surrounding fire.

FIRE FIGHTING PROCEDURES: Wear SCBA when fighting fires involving this as a precaution.

UNUSUAL FIRE AND EXPLOSION: None

HAZARDOUS COMBUSTION PRODUCTS: Carbon Monoxide, Carbon Dioxide, smoke, organic vapors from dried and heated

towels in fore conditions.

Section VI - Accidential Release Measures

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED:

Isolate damaged goods. Spills should be diluted with water, then absorbed with sand, clay, or earth. Dispose of saturated absorbent materials appropriately since spontaneous heating may occur.

Section VI - Handling and Storage

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE:

Isolate damaged goods. Spills should be diluted with water, then absorbed with sand, clay, or earth. Dispose of saturated absorbent materials appropriately since spontaneous heating may occur.

VENTILATION: Good general ventilation should be sufficient for most conditions.

RESPIRATORY PROTECTION: No respiratory protection should be needed.

PROTECTIVE CLOTHING: No precautions other than clean body covering clothing should be needed.

Gloves are recommended.

EYE PROTECTION: Use safety glasses or equivalent protection to avoid eye contact.

Section VIII - Exposure Controls / Personal Protection

HAZARDOUS COMPONENTS		NIOSH			ACGIH		0	SHA
Sodium Hypochlorite	TWA ppm	TWA mg/m3	STEL ppm	STEL mg/m3 10 IDLH	TLV/TWA ppm	TWA mg/m3	PEL ppm	PEL mg/m3

PERSONAL PROTECTION

RESPIRATORY PROTECTION: A respiratory protection program that meets OSHA's 29 CFR 1910.134, ANSI Z88.2 requirements

must be followed whenever workplace concentrations of this product in air exceed the published TLV's. Impervious rubber, or nitrile gloves are advised if prolonged or repeated exposure is likely.

EYE PROTECTION: Use safety glasses with protective side shields to avoid eye contact.

OTHER PROTECTIVE

PROTECTIVE GLOVES:

CLOTHING OR EQUIPMENT: Not required.

WORK / HYGIENIC PRACTICES: Wash with soap and water after contact.

microdot® Bleach Wipe

Section IX - Physical and Chemical Properties

APPEARANCE AND ODOR

Clear, colorless solution with a mild bleach odor.

Ha

10.0-12.0

BOILING POINT / BOILING RANGE

100"C (212°F)

FLASH POINT

> 1 01 QC (>215QF)

FLAMMABLE LIMITS

LFL:NF **UEL: NE**

AUTOIGNITION TEMPERATURE

VAPOR PRESSURE

17.5mmHg at 20 "C

SPECIFIC GRAVITY (WATER = 1)

1.015@25QC.

SOLUBILITY IN WATER

Completely soluble (for liquid).

VISCOSITY

NA

VAPOR DENSITY (AIR = 1)

NR

EVAPORATION RATE (WATER = 1)

Section IX - Stability and Reactivity

CONDITIONS TO AVOID: STABILITY

The product is stable under normal use and storage conditions.

Excessive heat and light exposure. Avoid contact with incompatible materials listed in the following section:

INCOMPATIBILITY (MATERIALS TO AVOID):

Avoid contact with acids, reducing agents, ammonia, or heavy metals such as nickel, cobalt, copper, and iron.

Reacts with other household chemicals such as toilet bowl cleaners, rust removers, vinegar, acids or ammonia containing products to produce hazardous gases, such as chlorine and other chlorinated species.

HAZARDOUS DECOMPOSITION OR BYPRO DUCTS:

Chlorine in the presence of acids.

HAZARDOUS POLYMERIZATION:

CONDITIONS TO AVOID:

Will not occur. None related to polymerization.

Section XI - Toxicological Information

HAZARDOUS COMPONENTS CAS# LD50 of Ingredient LC50 of Ingredient (Specify Species and Route) (Specify Species)

Sodium hypochlorite 7681-52-9 Oral, mouse: LD50 = 5800 mg/kg

In the case of medical emergency, contact your local poison control center.

Section XII - Ecological Information

No information is available for this blended product.

microdot® Bleach Wipe

Section XIII - Disposal Considerations

WASTE DISPOSAL METHOD:

Isolate damaged goods. Spills should be diluted with water, then absorbed with sand, clay, or earth. Dispose of saturated absorbent materials appropriately since spontaneous heating may occur. Landfill in a permitted waste disposal facility in accordance with all local, state, and federal regulations.

Section XIV - Transport Information

IATA HAZARD CLASS / Pack Group: Contact Manufacturer **PROPER SHIPPING NAME:** Not regulated

IMDG HAZARD CLASS: Contact Manufacturer **DOT HAZARD CLASS I Pack Group:** Not regulated

RID/ADR Dangerous Goods Code: Contact Manufacturer **REFERENCE:** 49CFR172, 49CFR173

UN TDG Class I Pack Group: Contact Manufacturer UN / NA IDENTIFICATION NUMBER: Not Regulated

Hazard Identification Number (HIN): Contact Manufacturer LABEL: Not Regulated

Refer to 49CFR172 and 49CFR173 for pertinent regulations **HAZARD SYMBOLS: NONE**

Note: Transportation information provided is for reference only. Client is urged to consult CFR 49 parts 100 - 177, IMDG, IATA, EU, United Nations TOG, and WHMIS (Canada) TOG information manuals for detailed regulations and exceptions covering specific container sizes, packaging materials and methods of shipping.

Section XV - Regulatory Information

TSCA (USA - Toxic Substance CONTRO ACT): Listed.

SARA TITLE III (USA - Superfund Amendments and Reauthorization Act):

Acute Health: Yes Chronic Health: Sudden Release of Pressure: No Fire: Nο

Reactive: Nο

313 REPORTABLE INGREDIENTS: None reportable.

ERCLA (USA - Comprehensive Response Compensation and Liability Act): CAS# 7681-52-9: 100 lb final RQ; 45.4

kg final RO. CAS# 1310-73-2: 1000 lb

final RQ; 454 kg final RQ.

Section XVI - Other Information

Legend: NA=Not Applicable	NE=Not established	NR=Not Reported	
HMIS HAZARD RATINGS	HEALTH: FLAMMABILITY: PHYSICAL HAZARD: PERSONAL PROTECTIVE: EQUIPMENT:	2 0 1 B Glasses and gloves	0 = INSIGNIFICANT 1 = SLIGHT 2 = MODERATE 3 = HIGH 4 = EXTREME
NFPA HAZARD RATINGS	HEALTH: FLAMMABILITY: REACTIVITY: SPECIAL HAZARDS::	0 0 1 OX	0 = NONE 1 = SLIGHT / LITTLE 2 = MODERATE 3 = HIGH / SERIOUS 4 = EXTREME

microdot® Minute Wipe

Product and Company Identification

NAIC 88494-1 (Wipe) Product identifier microdot® Minute Wipe Other means of identification

Disinfectant Recommended use None known Recommended restrictions

Manufacturer/Importer/Supplier/Distributor information

Manufacturer

Cambridge Sensors USA, LLC Company name

Address 10051 Bode Road

Plainfield, IL 60585 **United States**

1.877.374.4062 Telephone E-mail info@microdotcs.com

Emergency phone number 1.877.374.4062

Section II - Hazards Identification

Physical hazards Flammable liquids Category 2 **Health hazards** Acute toxicity, inhalation Category 4

Serious eye damage/eye irritation Category 2A Specific target organ toxicity, single exposure Category 3 respiratory tract irritation

Specific target organ toxicity, single exposure Category 3 narcotic effects

Environmental hazards Not classified. **OSHA** defined hazards Not classified.

Label elements



Danger Signal word

Highly flammable liquid and vapor. Causes serious eye irritation. Harmful if inhaled. May cause **Hazard statement**

respiratory irritation. May cause drowsiness or dizziness.

Precautionary statement

Prevention Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Keep container tightly

closed. Ground/bond container and receiving equipment. Use explosion-proof

electrical/ventilating/lighting equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Avoid breathing mist or vapor. Wash thoroughly after handling. Use only outdoors or in a well-ventilated area. Wear protective gloves/eye

protection/face protection.

Response In case of fire: Use appropriate media to extinguish. If on skin (or hair): Take off immediately all

contaminated clothing. Rinse skin with water/shower. If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. If inhaled: Remove person to fresh air and keep

comfortable for breathing. Call a poison center/doctor// if you feel unwell.

Storage Store in a well-ventilated place. Keep container tightly closed. Keep cool. Store locked up.

Disposal Dispose of contents/container in accordance with local/regional/national/international regulations.

Hazard(s) not otherwise None known.

classified (HNOC)

Supplemental information Although the product as a whole is in solid format, the product does not meet the OSHA HCS

definition of a flammable solid as per Appendix B to 1910.1200 - Physical Hazard Criteria, section B.7.1 and B. 7.2.

This is a registered EPA product. The product labeling is in compliance with EPA regulations and guidelines.

microdot® Minute Wipe

Section III - Composition/Information on Ingredients

Mixtures			
Chemical name	Common name and synonyms	CAS number	%
Ethanol		64-17-5	60.90

US GHS: The exact percentage (concentration) of composition has been withheld as a trade Composition comments secret in accordance with paragraph (i) of §1910.1200.

Section IV - First Aid Measures

If inhaled: Remove person to fresh air and keep comfortable for breathing. Call a poison Inhalation

center/doctor if you feel unwell.

If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. Skin contact

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and Eye contact easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Rinse mouth. Do not induce vomiting. Get medical attention if symptoms occur. Never give Ingestion

anything by mouth if victim is unconscious, or is convulsing. Obtain medical attention.

Most important Symptoms may include stinging, tearing, redness, swelling, and blurred vision. Vapors have a

narcotic effect and may cause headache, fatigue, dizziness and nausea. symptoms/effects, acute and

delayed

Indication of immediate Provide general supportive measures and treat symptomatically. In case of shortness of breath, medical attention and special give oxygen. Symptoms may be delayed. treatment needed

General information Ensure that medical personnel are aware of the material(s) involved, and take precautions to

protect themselves. Take off all contaminated clothing immediately. If you feel unwell, seek medical advice (show the label where possible). Show this safety data sheet to the doctor in attendance. Wash contaminated clothing before reuse. Keep away from sources of ignition. No smoking. Avoid contact with eyes and skin. Wear rubber gloves and safety glasses with side

shields. Keep out of reach of children.

Section V- Fire Fighting Measures

Suitable extinguishing media Dry chemical, CO2, water spray or regular foam.

Do not use water jet as an extinguisher, as this will spread the fire. Unsuitable extinguishing

media

Specific hazards arising from During fire, gases hazardous to health may be formed. the chemical

Special protective equipment Firefighters should wear full protective clothing including self contained breathing apparatus. and precautions for firefighters

Fire-fighting In case of fire and/or explosion do not breathe fumes. Move containers from fire area if you can do equipment/instructions so without risk.

Use standard firefighting procedures and consider the hazards of other involved materials. Specific methods

General fire hazards Solid containing flammable liquid

Section VI- Accidental Release Measures

Personal precautions, Remove all possible sources of ignition in the surrounding area. Wear appropriate protective equipment and clothing during clean-up. Do not breathe vapor. For personal protection, see protective equipment and section 8 of the SDS. emergency procedures

Methods and materials for Eliminate all ignition sources (no smoking, flares, sparks, or flames in immediate area). Take containment and cleaning up

precautionary measures against static discharge. Use only non-sparking tools. Pick up and discard towelette.

Avoid discharge into drains, water courses or onto the ground. Use appropriate containment to **Environmental precautions** avoid environmental contamination.

microdot® Minute Wipe

Section VII - Handling and Storage

Precautions for safe handling

Vapors may form explosive mixtures with air. Do not handle, store or open near an open flame, sources of heat or sources of ignition. Protect material from direct sunlight. Do not smoke. Minimize fire risks from flammable and combustible materials (including combustible dust and static accumulating liquids) or dangerous reactions with incompatible materials. Avoid breathing vapor. Avoid contact with eyes, skin and clothing. Avoid prolonged exposure. Use only outdoors or in a well-ventilated area. Wear appropriate personal protective equipment. Observe good industrial hygiene practices. Wash thoroughly after handling. Keep container tightly closed. When using do not eat or drink.

Conditions for safe storage, including any incompatibilities

Store locked up. Keep away from heat, sparks and open flame. Avoid spark promoters. Eliminate sources of ignition. Store in a closed container away from incompatible materials. Store in a cool, dry place out of direct sunlight. Store in a well-ventilated place. Store away from incompatible materials (see Section 10 of the SDS). Keep out of reach of children.

Section VIII - Exposure Controls/Personal Protection

Occupational exposure limits

US. OSHA Table Z-1 Limits for Air	Contaminants (29 CFR 1910.1000)
Componento	Type

Components	туре	value	
Ethanol (CAS 64-17-5)	PEL	1900 mg/m3	_
		1000 ppm	
US. ACGIH Threshold Limit Valu	ies		
Components	Type	Value	
Ethanol (CAS 64-17-5)	STEL	1000 ppm	
US. NIOSH: Pocket Guide to Che	emical Hazards		
Components	Туре	Value	
Ethanol (CAS 64-17-5)	TWA	1900 mg/m3	
		1000 ppm	

Biological limit values

No biological exposure limits noted for the ingredient(s).

controls

Appropriate engineering

Good general ventilation (typically 10 air changes per hour) should be used. Ventilation rates should be matched to conditions. If applicable, use process enclosures, local exhaust ventilation, or other engineering controls to maintain airborne levels below recommended exposure limits. If exposure limits have not been established, maintain airborne levels to an acceptable level.

Individual protection measures, such as personal protective equipment

Eye/face protection Skin protection

Wear safety glasses with side shields.

Hand protection

Rubber gloves. Confirm with a reputable supplier first.

Other

Wear suitable protective clothing. Wear appropriate chemical resistant clothing. As required by

employer code.

Respiratory protection

Where exposure guideline levels may be exceeded, use an approved NIOSH respirator.

Thermal hazards

Not applicable.

General hygiene considerations

When using do not smoke. Always observe good personal hygiene measures, such as washing after handling the material and before eating, drinking, and/or smoking. Routinely wash work

clothing and protective equipment to remove contaminants.

microdot® Minute Wipe

Section IX - Physical and Chemical Properties

Appearance Clear Liquid saturated on wipe

Physical state Liquid.

Form Liquid saturated on wipe

ColorColorlessOdorAlcoholOdor thresholdNot available.pH3.7 @ 77°F (liquid)Melting point/freezing pointNot available.

range

Pour pointNot available.Specific gravityNot available.Partition coefficientNot available.

(n-octanol/water)

Initial boiling point and boiling

Flash point 65.3 °F (18.5 °C) Tag Closed Cup (liquid)

Not available.

Evaporation rate Not available.
Flammability (solid, gas) Not applicable.
Upper/lower flammability or explosive limits

Flammability limit - lower

Not available.

(%)

Flammability limit - upper

Not available.

(%)

Explosive limit - lower (%) Not available.

Explosive limit - upper (%) Not available.

Vapor pressure Not available.

Vapor density Not available.

Relative density Not available.

Solubility(ies) Not available.

Auto-ignition temperature Not available.

Viscosity 3.23 cSt @ 68°F (liquid)

Other information

Decomposition temperature

Density 0.86 (liquid)

Section X - Stability and Reactivity

Reactivity This product may react with strong oxidizing agents.

Not available.

Possibility of hazardous No dangerous reaction known under conditions of normal use.

reactions

Chemical stability Stable under recommended storage conditions.

Conditions to avoid Avoid heat, sparks, open flames and other ignition sources. Avoid temperatures exceeding the

flash point. Do not mix with other chemicals.

Incompatible materials Strong oxidizing agents.

Hazardous decomposition

products

May include and are not limited to: Oxides of carbon.

microdot® Minute Wipe

Section XI - Toxicological Information

Information on likely routes of exposure

Ingestion Expected to be a low ingestion hazard.

Vapors have a narcotic effect and may cause headache, fatigue, dizziness and nausea. Harmful if inhaled. Prolonged inhalation may be harmful. Inhalation

Skin contact No adverse effects due to skin contact are expected. Non-irritating based on test data.

Eye contact Causes serious eye irritation.

Symptoms may include stinging, tearing, redness, swelling, and blurred vision. Symptoms of Symptoms related to the overexposure may be headache, dizziness, tiredness, nausea and vomiting.

physical, chemical and toxicological characteristics

Information on toxicological effects

Harmful if inhaled. Narcotic effects. May cause respiratory irritation. Acute toxicity **Test Results**

Product Species NAIC 88494-1 (Wipe) (CAS Mixture)

Acute

Dermal

LD50 Rat

Rat

> 5000 mg/kg, Tested

Inhalation LC50

> 2.1 mg/l, Tested

Oral LD50

Rat

> 5000 mg/kg, Tested

Components Species Test Results

Ethanol (CAS 64-17-5)

Acute

Dermal

Rabbit > 15800 mg/kg

Inhalation

LC50 Mouse

Rat

39 mg/l, 4 Hours 31623 ppm, 4 Hours

3450 ma/ka

20000 ppm, 10 Hours

Oral LD50

Dog 5500 mg/kg Guinea pig 5600 mg/kg

Mouse Rat

7060 mg/kg Skin corrosion/irritation Prolonged skin contact may cause temporary irritation.

Exposure minutes Not available. Erythema value Not available.

Oedema value Not available.

Serious eye damage/eye

Causes serious eye irritation.

Corneal opacity value Not available. Iris lesion value Not available Conjunctival reddening Not available.

value

Conjunctival oedema value Not available. Recover days Not available.

Respiratory or skin sensitization

Respiratory sensitization

Skin sensitization This product is not expected to cause skin sensitization.

Germ cell mutagenicity Non-hazardous by OSHA criteria. Carcinogenicity Non-hazardous by OSHA criteria.

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Not listed.

Reproductive toxicity Non-hazardous by OSHA criteria. Respiratory tract irritation. Narcotic effects.

Specific target organ toxicity single exposure

Not classified.

Specific target organ toxicity repeated exposure

Aspiration hazard Not available.

Prolonged inhalation may be harmful. Chronic effects

Not available. **Further information**

microdot® Minute Wipe

Section XII - Ecological Information

Eco	otoxicity	See below		
	Components		Species	Test Results
	Ethanol (CAS 64-17-5)			
	Crustacea	EC50	Daphnia	11744.5 mg/L, 48 Hours
	Aquatic			
	Crustacea	EC50	Water flea (Daphnia magna)	7.7 - 11.2 mg/l, 48 hours
	Fish	LC50	Fathead minnow (Pimephales promelas)	> 100 mg/l, 96 hours

No data is available on the degradability of this product. Persistence and degradability

Bioaccumulative potential No data available. Partition coefficient n-octanol / water (log Kow)

Ethanol -0.31

Mobility in soil No data available. Not available. Mobility in general

Other adverse effects No other adverse environmental effects (e.g. ozone depletion, photochemical ozone creation potential, endocrine disruption, global warming potential) are expected from this component.

Section XIII - Disposal Considerations

Disposal instructions Collect and reclaim or dispose in sealed containers at licensed waste disposal site. This material and its container must be disposed of as hazardous waste. Incinerate the material under controlled conditions in an approved incinerator. Do not incinerate sealed containers. Do not allow this material to drain into sewers/water supplies. Do not contaminate ponds, waterways or ditches with

chemical or used container. If discarded, this product is considered a RCRA ignitable waste, D001. Dispose of contents/container in accordance with local/regional/national/international regulations.

Local disposal regulations Dispose in accordance with all applicable regulations. Hazardous waste code

D001: Waste Flammable material with a flash point <140 F

The waste code should be assigned in discussion between the user, the producer and the waste

disposal company.

Waste from residues / unused

products

Dispose of in accordance with local regulations. Empty containers or liners may retain some product residues. This material and its container must be disposed of in a safe manner (see:

Disposal instructions).

Empty containers should be taken to an approved waste handling site for recycling or disposal. Contaminated packaging

Since emptied containers may retain product residue, follow label warnings even after container is

Section XIV - Transport Information

U.S. Department of Transportation (DOT)

Basic shipping requirements:

UN number

Proper shipping name Flammable liquids, toxic, n.o.s. (Ethanol)

Hazard class Limited Quantity - US

Packing group Ш

IB2, T7, TP2, TP13 Special provisions

< 0.3 Gal-Consumer Commodity ORM-D/Limited Quantity Packaging exceptions

202 Packaging non bulk Packaging bulk 243

DOT



microdot® Minute Wipe

Section XII - Regulatory Information

US federal regulations This product is a "Hazardous Chemical" as defined by the OSHA Hazard Communication

Standard, 29 CFR 1910.1200.

TSCA Section 12(b) Export Notification (40 CFR 707, Subpt. D)

Not regulated.

CERCLA Hazardous Substance List (40 CFR 302.4)

US. OSHA Specifically Regulated Substances (29 CFR 1910.1001-1050)

Superfund Amendments and Reauthorization Act of 1986 (SARA)

Immediate Hazard - Yes **Hazard categories**

Delayed Hazard - No Fire Hazard - Yes Pressure Hazard - No Reactivity Hazard - No

SARA 302 Extremely

hazardous substance

SARA 311/312 Hazardous

chemical

SARA 313 (TRI reporting)

Not regulated.

Other federal regulations

Clean Air Act (CAA) Section 112 Hazardous Air Pollutants (HAPs) List

Clean Air Act (CAA) Section 112(r) Accidental Release Prevention (40 CFR 68.130)

Not regulated.

Safe Drinking Water Act

Not regulated.

(SDWA)

Food and Drug Not regulated.

Administration (FDA)

Ethanol is not considered to be a Prop 65 component for this product. The product is not an US state regulations

alcoholic beverage for consumption.

US - California Hazardous Substances (Director's): Listed substance

Ethanol (CAS 64-17-5) Listed.

US - California Proposition 65 - Carcinogens & Reproductive Toxicity (CRT): Listed substance

Ethanol (CAS 64-17-5) Listed.

US - Illinois Chemical Safety Act: Listed substance US - Louisiana Spill Reporting: Listed substance

Ethanol (CAS 64-17-5) Listed.

Ethanol (CAS 64-17-5) Listed.

US - Minnesota Haz Subs: Listed substance

Ethanol (CAS 64-17-5) Listed. US - New Jersey RTK - Substances: Listed substance

Ethanol (CAS 64-17-5) Listed. US - Texas Effects Screening Levels: Listed substance

Ethanol (CAS 64-17-5) Listed.

US. Massachusetts RTK - Substance List

Ethanol (CAS 64-17-5) Listed

US. Pennsylvania RTK - Hazardous Substances

Ethanol (CAS 64-17-5) Listed.

US. Rhode Island RTK

Not regulated.

On inventory (yes/no)* Country(s) or region Inventory name

United States & Puerto Rico Toxic Substances Control Act (TSCA) Inventory

*A "Yes" indicates this product complies with the inventory requirements administered by the governing country(s)

Yes

microdot® Minute Wipe

Section XII - Other Information


LEGEND	
Severe	4
Serious	3
Moderate	2
Slight	1
Minimal	0

HEALTH * 1

FLAMMABILITY 3

PHYSICAL HAZARD 0

PERSONAL X



Disclaimer

The information in the sheet was written based on the best knowledge and experience currently available. Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since conditions of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirements of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document.

Issue date 30-June-2015

Further information For an updated SDS, please contact the supplier/manufacturer listed on the first page of the

document

Other information This Safety Data Sheet was prepared to comply with the current OSHA Hazard Communication

Standard (HCS) adoption of the Globally Harmonized System of Classification and Labeling of

Chemicals (GHS).

Prepared by Cambridge Sensors USA, LLC Phone: 1.877.374.4062

TROUBLESHOOTING THE microdot® Xtra METER



The blood glucose level is higher than 525 mg/dL.

This message indicates very high blood sugar. You should recheck the blood glucose level. If "Hi" again, follow facility protocol.

™ 8:08 8-8

The blood glucose level is lower than 20 mg/dL.

This message indicates very low blood sugar. You should recheck the blood glucose level. If "Lo" again, follow facility protocol.

E-1 10 7

Error message that indicates that there is a problem with the meter, e.g., measurement error (time out, overflow, offset) or temperature out of range.

Review the instructions and try again with a new test strip. If the problem persists, contact Customer Service.

201

Error message that indicates that there may be a problem with the test strip, e.g, the test strip may be damaged, moved, or removed during testing, or inserted improperly.

Check the test strip for damage and retest as necessary. Repeat the test. If the error message appears again, contact Customer Service.

E-3 311

Error message could be caused by a used or damaged test strip.

Repeat the test with a new strip. If the error message appears again, contact Customer Service.

E-5 50 3

Error message indicates serial communications error.

If error persists, contact Customer Service.

The battery sign appears on the display with the unit of measurement. The power of the battery is getting low. You can complete about 50 more tests from the time this symbol first appears.

Test results will be accurate but replace the battery as soon as possible. Battery type is CR2032 3 volt.

microdot xtra

Cleaning & Disinfecting the microdot® Xtra Meter

Prior to cleaning & disinfecting remove the orange, protective microdot® Meter Sleeve. Cambridge Sensors USA, LLC recommends cleaning & disinfecting the microdot® Xtra Meter using an EPA Registered 1:10 dilution of .525% - .615% sodium hypochlorite, or an EPA-registered tuberculocidal disinfectant wipe. The recommendation is in accordance with C.D.C. Guidelines for Disinfection and Sterilization.

Cambridge Sensors USA, LLC offers two ready to use products to meet cleaning & disinfection requirements. The ready to use microdot[®] Bleach Wipe and microdot[®] Minute Wipe. The microdot[®] Bleach Wipe is premoistened with a 1:10 dilution. The microdot® Minute Wipe is an EPA-registered tuberculocidal.

For easy cleaning & disinfecting Cambridge Sensors USA, LLC offers the microdot® Disinfection System which includes the user friendly microdot® Meter Case & Timer to insure proper dwell time.

Note: Over exposure (contact time) may damage meter. When cleaning & disinfecting the microdot® Xtra Meter take extreme care not to get liquid in the test strip dock or key code parts.

Cambridge Sensors USA, LLC

Customer Service available from from 8am - 5pm CST Monday - Friday Toll free 1.877.374.4062 | www.microdotcs.com

Non-Bleach Wipe Protocol to Disinfect your microdot® Xtra Glucometer

Proper Guidelines for Disinfection of microdot® Xtra Glucometer with Disinfection Case and Timer

Read the microdot® Minute Wipe (EPA Reg No. 88494-2-88459) label and follow directions for use. Wipe is for the exterior surfaces of blood glucose meters and is not for the use on the needle or the monitor.

- Always use personal protective equipment as specified on the microdot[®] Minute Wipe label.
- Thoroughly clean gross filth and heavy soil from surface of microdot® Xtra Glucometer to be disinfected.



- 1. Remove a pre-saturated 6" x 6.75" microdot® Minute Wipe.
- 2. Thoroughly wipe the microdot® Xtra Glucometer surface to be disinfected.
- 3. Wrap the Glucometer with the microdot® Minute Wipe.



4. Place the wrapped microdot® Xtra Glucometer face down inside the microdot® Disinfection Case.



- 5. Close disinfection case lid and activate 1 minute timer.
- 6. Allow the microdot® Xtra Glucometer to remain in contact with the microdot® Minute Wipe for 1 minute.
- 7. Dispose of wipe in trash after use. Do not flush wipe in the toilet.
- 9. Dispose of the non-refillable empty canister according to state and local authorities guidelines as allowed by the microdot® Minute Wipe label.

NOTE: Over exposure (contact time) may damage microdot® Xtra Glucometer.

Bleach Wipe Protocol to Disinfect your microdot® Xtra Glucometer

Proper Guidelines for Disinfection of microdot® Xtra Glucometer with disinfection case and timer

Read the microdot® Bleach Wipe (EPA Reg. No. 37549-1-88459) label and follow directions for use. microdot® Bleach Wipe is for the exterior surfaces of blood glucose meters and is not for the use on the needle or the monitor.

- Always use personal protective equipment as specified on the microdot® Bleach Wipe label.
- Thoroughly clean gross filth and heavy soil from surface of microdot® Xtra Glucometer to be disinfected.
- 1. Open microdot® Bleach Wipe pop-up canister. The wipes are pre-saturated with a sodium hypochlorite (bleach) hospital-use solution.



- 2. Remove a pre-saturated 6" x 6" wipe.
- 3. Thoroughly wipe the microdot® Xtra Glucometer surface to be disinfected.
- 4. Wrap the Glucometer with the microdot® Bleach Wipe.



5. Place the wrapped microdot® Xtra Glucometer face down inside the microdot® Disinfection Case.



- 6. Close disinfection case lid and activate 3 minute timer.
- 7. Allow the microdot® Xtra Glucometer to remain in contact with the bleach wipe for 3 minutes.
- 8. Dispose of wipe in trash after use. Do not flush wipe in the toilet.
- 9. Dispose of the non-refillable empty canister according to state and local authorities guidelines as allowed by the microdot® Bleach Wipe label.

NOTE: Over exposure (contact time) may damage microdot® Xtra Glucometer.



CS° CAMBRIDGE SENSORS USA, LLC

CAMBRIDGE SENSORS USA, LLC 10051 Bode Rd. Plainfield, IL 60585

Customer Service: toll free: 877.374.4062 Monday-Friday 9am-5pm CST

www.microdotcs.com