

CAT[®] TOURNIQUET *Combat Application Tourniquet*

REV030218

Proven to be 100% Effective...

The patented C-A-T[®] is a true one-handed tourniquet proven to be 100% effective by the U.S. Army's Institute of Surgical Research. Tests proved that the C-A-T[®] completely occluded blood flow of an extremity in the event of a traumatic wound with significant hemorrhage.

The C-A-T[®] utilizes a durable windlass system with a patented free-moving internal band providing true circumferential pressure to the extremity. Once adequately tightened, bleeding will cease and the windlass is locked into place. A hook and loop windlass retention strap is then applied, securing the windlass to maintain pressure during casualty evacuation. The C-A-T[®]'s unique dual securing system avoids the use of screws and clips which can become difficult to operate under survival stress or where fine motor skills are compromised.

The patented Combat Application Tourniquet[®] (C-A-T[®]) featuring the proprietary red tip design and the mechanical advantage of a band within a band has been the Official Tourniquet of the U.S. Army since 2005. The innovations in the C-A-T[®] are protected by U.S. Patent Nos. 7,842,067 and 7,892,253.

Features:

- Official Tourniquet of the U.S. Army
- Proven to be 100% effective in occluding blood flow in both upper & lower extremities by the U.S. Army's Institute of Surgical Research
- NEW Single Routing Buckle for faster application, decreased blood loss, effective slack removal, fewer windlass turns, and simplified training with single protocol application standards
- Includes reinforced windlass, stabilization plate, windlass clip and writeable windlass security strap
- Featuring NAR's signature Red Tip Technology[®] elliptical tab providing added visual clues during application

Dimensions:

- Packaged: L 6.5 in. x W 2.4 in. x D 1.5 in.
- Open Length: 37.5 in.
- Weight: 2.7 oz

*The official
tourniquet of the
U.S. Army*



C-A-T[®] Tourniquet Gen 7

Color	Item #	NSN#
BLK	30-0001	6515-01-521-7976
ORG	30-0023	6515-01-612-7129
BLUE-Trainer	30-0033	6910-01-560-2972



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION[®]**