



Nitrile

An effective barrier against blood-borne pathogens

- Advanced non-latex nitrile polymer offers superior stretch and comfort
- Meets or exceeds the emergency medical glove requirements of NFPA 1999
- Fentanyl Permeation Testing: zero breakthrough with Fentanyl Citrate 100mcg/2ml for 240 minutes
- Non-stick exterior surface prevents tapes from sticking to the glove
- AQL of 0.65 exceeding industry standards
- Human skin pH 5.5 balanced coating for reduced risk of irritation and perspiration



Industries

- Emergency Medical Services
- Healthcare
- Law Enforcement

Recommended For

- Administering drugs
- Standard patient examinations
- Emergency services



TECHNICAL DATA SHEET

Product Information	
External Glove Surface	Fully Textured
Audit Standards	ISO 13485, ISO 9001
Product Certification	Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country Of Origin	China
Product Segmentation	High Risk
Antistatic	Not Tested


PHYSICAL PROPERTIES

	Typical Values		Testing Method
	BEFORE AGING	AFTER AGING	
Length (mm/inches)	300 / 11.8		ASTM D3578,ASTM D3577
Freedom from Holes (Inspection level I)	0.65 AQL		ASTM D5151
Palm Thickness (mm/mils)	0.16 / 6.3		ASTM D3578,ASTM D3577
Finger Thickness (mm/mils)	0.18 / 7.1		ASTM D3578,ASTM D3577
Ultimate Tensile Strength (MPa)	20	18	ASTM D412 & D573
Elongation at Break (%)	550	500	ASTM D412 & D573

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	XXL (10.5 - 11)	XXXL (11.5 - 12)
Product Code	FL12-1	FL12-2	FL12-3	FL12-4	FL12-5	FL12-6

Performance Standards and Regulatory Compliance



ISO 13485 ISO 9001

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ansell PYC
Ten. +7 495 258 13 16

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking © 2022 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Life-Assist
Helping Heroes Save Lives

Ansell