

Nitrile Powder-Free Exam Gloves

Attributable to the continuous advancement of conventional nitrile gloves, the Amsino Nitrile Glove Series provide exceptional barrier protection while enhancing comfort and dexterity during prolonged procedures.

Unique Properties

- Excellent film strength to minimize failures
- Excellent puncture resistance
- Good conformability reducing hand fatigue
- Textured fingertips provide good grip strength
- Ease in donning

Allergy Care

Does not contain natural rubber latex protein

Product Attributes

- FDA 510(k) Medical Use
- Textured Fingertips
- Beaded Cuff 9.5 inches
- Barrier Integrity: AQL 1.5
- Residual Powder Content: < 2.0 mg/glove
- ASTM D6124



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768

800.632.6746 • fax: 909.626.3888 • info@amsino.com



Product Specification

(Meets ASTM D6319 and EN 455)

Dimension (mm)	Spec	Typical		Physical Properties	
		mm	mil		
Palm	Min. 0.05	0.06	2.3	Tensile Strength (MPa)	Min. 14
Finger	Min. 0.05	0.09	3.5	Elongation (%)	Min. 500
Length	Min. 230	235		Force at Break (N)	Min. 6.0
Color	Blue			Tensile Strength (MPa)	Min. 14
				Elongation (%)	Min. 400
				Force at Break (N)	Min. 6.0

Item Code	Size	Description	Quantity/Unit	Case Dimension
ASCS38S	Small	Nitrile Powder Free Examination Glove, Blue	100 Pieces/Box 10 Boxes/Case 1000 Pieces/Case	315 x 250 x 250 mm
ASCS38M	Medium			Box Dimension
ASCS38L	Large			
ASCS38XL	X-Large			

Call Amsino sales reps for pricing, lead time and any other order information.

To view our complete product line, visit us at www.amsino.com

Disclaimer: It is important to review the chemical test report(s) to determine the required level of protection as result may vary for each glove. Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2020, Amsino International, Inc. All rights reserved.