

# Ansell Fentanyl Tested Products -MICROFLEX

Brand	Product Name	Style No.	Testing For Fentanyl & Gastric Acid <small>(min. &gt;240 minutes)</small>
MICROFLEX	LifeStar EC	LSE-104	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	Supreno EC	SEC-375	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM <b>D6978-05 pass/fail mark is &gt;240 minutes</b> ) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	Supreno SE	SU-690	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)
MICROFLEX	MidKnight XTRA	93-862	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	MidKnight TOUCH	93-732	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	Blaze	N48	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)
MICROFLEX	Freeform SE Freeform EC	FFS-700 FFE-775	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)
MICROFLEX	MidKnight	MK-296	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	Cobalt	N19	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)
Micro-Touch	Nitrafree Blue		<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes)
Micro-Touch	Royal Blue	313029	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol
MICROFLEX	UltraSense	US-220	<b>Static Testing</b> used for both Fentanyl & Gastric Acid (ASTM D6978-05 pass/fail mark is >240 minutes) <b>ACPP (Dynamic Testing Protocol)*</b> used for Fentanyl testing - Results showed no permeation after 30 min. of dynamic testing protocol