

# BURNFREE PRODUCTS

## 12" x 12" (30cm x 30cm) Sterile Burn & Wound Dressing



- Provides immediate pain relief by covering exposed nerve endings
- Rapidly cools and soothes the burn
- Draws the heat out and stops the burn progression
- Reduces swelling, redness, and additional skin and tissue damage
- Medical grade ingredients
- Portable, safe, and easy to use
- Sterile dressing

<b>INTL/US: CASE DIMENSIONS</b>	QUANTITY: 20 dressings WEIGHT: 12.2 pounds LENGTH: 16 inches WIDTH: 8 inches HEIGHT: 8 inches
<b>PRODUCT CODES</b>	UPC: 702135412122 INTL/EU: I-1212-20 INTL/US: 1212-20
<b>U.S. PATENT</b>	5,958,420
<b>INGREDIENTS</b>	Purified Water, Melaleuca Alternifolia, Emulsifiers and Preservatives

BurnFree solves the immediate first aid needs for all types of burns by relieving the pain, rapidly removing the heat, and cooling the burn.

BurnFree Sterile Dressings are designed for first aid use on all types of burn injuries. The dressings are constructed of a special, soft, open-cell foam, and are elastic enough to accommodate possible swelling or the need to flex joints. BurnFree dressings are non-adherent and the inert foam will not shed threads, strands, or fibers that could become embedded in the wound. The pliable foam holds 50 times its weight in gel. This means that more BurnFree gel will stay in contact with the burn.

The 12" x 12" (30cm x 30cm) BurnFree Sterile Dressing can be used for larger burns on arms, legs, shoulders, back, and face. For burns on the face, the foam can easily be cut or torn to accommodate for eyes, nose, or mouth.



**EXURO**<sup>LLC</sup>  
medical<sup>®</sup>

**BURNFREE** A Division of Exuro Medical LLC  
310 East 300 South, Salt Lake City, UT 84111  
Ph: 801-569-9090 | Fax: 801-569-3733  
info@exuromedical.com | exuromedical.com