

Casualty Care Rescue Randy product profile

Product Name	Casualty Care Rescue Randy - powered by Strategic Operations Hyper-Realistic® technology
Addresses Training Needs/Curricula	TCCC - Tactical Combat Casualty Care for all service members TECC - Tactical Emergency Casualty Care for all EMS/EMT/LEO professionals
Configuration	<ul style="list-style-type: none"> • Full-body manikin with articulated joints • Weight-distributed to mimic the human anatomy • Holds 3-4 liters of blood • Pumping system with regulator simulating human blood flow pressure (2-3) psi • Manual foot operated blood pumping system allows 3 wounds to bleed simultaneously • Optional Electric Deluxe blood pumping system allowing 4 wounds to bleed simultaneously • Comes with blood powder (makes 5 gallons) • Three (3) repairable neck skins • Ten (10) repairable tracheas • One (1) repair kit
Dimensions	Male at 72" tall x 25" wide x 12" deep (165.1 cm x 63.5 cm x 30.48 cm)
Product Ship Weight	80 lbs. Product Weight
Product Code/SKU	149-6075



Nasco Healthcare Inc.
16 Simulaids Drive
Saugerties, NY, 12477 USA
Phone: 1.800.431.4310
Email: info@nascohealthcare.com
www.nascohealthcare.com



Casualty Care Rescue Randy



Ready to treat. Ready to act.

Introducing Casualty Care Rescue Randy


Full-body manikin that is uniquely designed for enhanced realistic training on the procedures that treat the 3 most preventable causes of death:

1. Massive Bleeding
2. Tension Pneumothorax
3. Airway Obstruction

- Compliant with TCCC/TECC scenarios
- Rugged, reliable, and hyper-realistic
- 100% hand-made in the USA



Casualty Care Rescue Randy:
Life-like training solutions for
real-life situations.



This manikin is intended to be used as an instructional aide in the teaching and demonstrating of live rescue training. Actual product may vary slightly from photo. Nasco reserves the right to change product color, materials, supplies, or function as needed.