

SHOULDER HARNESS RESTRAINT SYSTEMS

The **Shoulder Harness Restraint System** is made up of two adjustable shoulder straps, each 60" long. And one 5' double adjusting torso strap. Metal Push Button buckle attaches into L-Links to secure on patient.

Item #11170C, 31170C, 51670C & 61170C fits Stryker Cot 6500 Series

- Available in **Nylon, Impervious, BioThane® G1** and **BioThane® G2** Webbing.

NYLON

Item # 11160

Nylon Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5½' Torso Strap**

Item # 11153

Nylon Shoulder Harness Restraint System, Torso Strap only, each, **5½' Torso Strap**

Item # 11154

Nylon Shoulder Harness Restraint System, Shoulder Straps only, pair, **4½' Shoulder Straps**

Item #11170

Nylon Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7½' Torso Strap**

Item #11170C

Nylon Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

IMPERVIOUS

Item # 31160

Impervious Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5½' Torso Strap**

Item # 31153

Impervious Shoulder Harness Restraint System, Torso Strap only, each, **5½' Torso Strap**

Item # 31154

Impervious Shoulder Harness Restraint System, Shoulder Straps only, pair, **4½' Shoulder Straps**

Item #31170

Impervious Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7½' Torso Strap**

Item #31170C

Impervious Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**



BIOTHANE® G1

Item # 51660

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5½' Torso Strap**

Item # 51653

BioThane® G1 Shoulder Harness Restraint System, Torso Strap only, each, **5½' Torso Strap**

Item # 51654

BioThane® G1 Shoulder Harness Restraint System, Shoulder Straps only, pair, **4½' Shoulder Straps**

Item #51670

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7½' Torso Strap**

Item #51670C

BioThane® G1 Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

BIOTHANE® G2

Item # 61160

BioThane® G2 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 5½' Torso Strap**

Item # 61153

BioThane® G2 Shoulder Harness Restraint System, Torso Strap only, each, **5½' Torso Strap**

Item # 61154

BioThane® G2 Shoulder Harness Restraint System, Shoulder Straps only, pair, **4½' Shoulder Straps**

Item #61170

BioThane® G2 Shoulder Harness Restraint System, Complete Set, **4½' Shoulder Straps, 7½' Torso Strap**

Item #61170C

BioThane® G2 Shoulder Harness Restraint System, Complete Set, **6½' Shoulder Straps, 6' Torso Strap**

SHOULDER HARNESS RESTRAINT WITH FIXED SHOULDER HOLDER



The **Fixed Shoulder Holder** is made with an 8" X 2" nylon strap sewn to both adjustable shoulder straps to prevent straps from falling off patients shoulders.

- Available in **Nylon, Impervious, BioThane® G1** and **BioThane® G2** Webbing.

Item # 11160F

Nylon Shoulder Harness System with Fixed Shoulder Holder, Complete Set

Item # 31160F

Impervious Shoulder Harness System with Fixed Shoulder Holder, Complete Set

Item # 51660F

BIOTHANE® G1 Shoulder Harness System with Fixed Shoulder Holder, Complete Set

Item # 61160F

BIOTHANE® G2 Shoulder Harness System with Fixed Shoulder Holder, Complete Set

