



Yankauer Handles with Suction Tubing


Dynarex dependable Yankauer handles with bulb and straight tips and suction tubing

- For convenient, affordable respiratory suctioning procedures
- Individually sealed in easy-open, peel-strip packages
- Medical grade PVC
- Sterile suction tool for removing particulates and viscous material
- Bulb tip to minimize tissue trauma
- Ribbed connector for secure attachment to suction tubing
- Blue-tinted, transparent material with slip-resistant handle


Yankauer Suction Handle
Straight Tip/Non-Vented
Reorder No. 4695

Yankauer Suction Handle
Bulb Tip/Non-Vented
Reorder No. 4691

Bulb Tip



Yankauer with Suction
Tubing – Sterile
Reorder No. 4698



Vent

Yankauer Suction Handle
Straight Tip/Vented
Reorder No. 4694

Yankauer Suction Handle
Bulb Tip/Vented
Reorder No. 4690

Yankauer Suction Tubing
Non-Conductive – Sterile
w/Straw Connectors
Reorder No. 4693

Reorder No.	Description	Quantity
4690	Yankauer Suction Handle – Bulb Tip/Vented	50/Cs
4691	Yankauer Suction Handle – Bulb Tip/Non-Vented	50/Cs
4693	Suction Tubing, Non Conductive – Sterile w/Straw Connector 9/32" x 6ft	50/Cs
4694	Yankauer Suction Handle – Straight Tip/Vented	50/Cs
4695	Yankauer Suction Handle – Straight Tip/Non-Vented	50/Cs
4698	Yankauer with Suction Tubing – Bulb Tip/Non-Vented 1/4" x 6ft	50/Cs

* Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.



Also see our website
for our Suction Canisters
(Reorder No. 4670 - 4672)

For more information contact your local Dynarex® distributor sales representative.
Or call **888-396-2739** to find your nearest distributor.

Already mobile?
Scan your smart phone here for more information.



Find us on:

