



Products / Traditional Wound Care / Burn Products / Burn Dressing / Burn Gauze Pad 18" x 18" 20" x 10" 1-Ply



Burn Gauze Pad 18" x 18" 20" x 10" 1-Ply

9170

100/CS

Product Features

100% Woven cotton gauze with thick ply

Provides superior absorption and a high level of padding and protection to site

Ideal for debridement of wounds

Mesh cut

Attributes

Brand Name	Dukal
Sterility	Sterile
Size	18" x 18" 20" x 10"
Ply	1-Ply
Contains Latex	Not made with natural rubber latex
HCPC Code	A6404
Weight	10.95

Packaging

Packaging	100/CS
Case Packaging	100/CS
Total Case Quantity	100
Item Packaging Barcode	(01)10665973001170
Case Packaging Barcode	(01)20665973001177

Dukal Corporation
2 Fleetwood Ct.
Ronkonkoma, NY 11779

Phone: 1 (631) 656-3800
Toll Free: 1 (800) 243-0741
Fax: 1 (631) 656-3810

