

## TECHNICAL DATA SHEET

<b>Graham Medical Item #:</b>	86792
<b>Brand:</b>	Graham Medical
<b>Product Type:</b>	Protective Gown
<b>Country of Origin:</b>	China
<b>Harmonized Code:</b>	6210.10.5000
<b>Color:</b>	Blue
<b>Material:</b>	LDPE Film
<b>Protection Level:</b>	AAMI Level 2
<b>Chemotherapy Tested:</b>	No
<b>Sterility:</b>	NonSterile
<b>User:</b>	Adult
<b>Usage:</b>	Disposable
<b>Back Style:</b>	Open Back
<b>Sleeve Style:</b>	Long Sleeves
<b>Cuff Style:</b>	Thumb-Loop
<b>Neck Closure:</b>	Over-The-Head, Easy-Tear Removal
<b>Fastening Type:</b>	Integrated Waist Tie
<b>Gender:</b>	Unisex
<b>Size:</b>	42"W x 46"L - One Size Fits Most
<b>Latex-Free:</b>	Yes
<b>UNSPSC Code:</b>	42131612
<b>Packaging:</b>	One Gown Per Bag



Revised: 07/10/20