

**PURPLE NITRILE\***

Exam Gloves



**Life-Assist**  
Helping Heroes Save Lives

# MAXIMUM PROTECTION FOR ALL RISKS

## THE WORRY-FREE CHOICE...DESIGNED FOR USE BY EVERYONE, EVERYWHERE

HALYARD\* PURPLE NITRILE\* Exam Gloves are often the glove of choice for healthcare workers who need proven barrier protection from bacteria, blood-borne viruses and chemicals during procedures where risk of fluid exposure is moderate to high. PURPLE NITRILE\* gloves offer the highest level of chemical protection available in Halyard's line of nitrile exam gloves. They are especially appropriate in dealing with unknown risks.

These durable gloves are available in 9.5 and 12 inch lengths, sterile and non-sterile versions, to provide a full range of protection options. 12-inch PURPLE NITRILE-XTRA\* Gloves are NFPA 1999-2008 Certified (National Fire Protection Association).

### TRUSTED PERFORMANCE

Healthcare professionals in the ER, ICU, Lab, and Oncology Units of acute care facilities, as well as EMS staff, Dental clinicians, and Cancer Center specialists all rely on PURPLE NITRILE\* brand of Exam Gloves. When they see that unique purple color, clinicians know they have protection they can trust, and that gives them the confidence to provide the best care to every patient.

### TESTED AGAINST 29 CHEMOTHERAPY DRUGS

PURPLE NITRILE\* Exam Gloves are an excellent choice for protection against chemical exposure. They have been tested against 29 chemotherapy drugs per the ASTM D6978-05 "Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs," as well as common chemicals found in the healthcare environment. Please ask your Halyard Representative for a complete list of the various levels of protection that PURPLE NITRILE\* exam gloves provide.



## Product Specifications

Gauge Thickness Measurements	MM	MIL
Middle Finger:	.15	5.9
Palm:	.12	4.7
Cuff:	.09	3.5
Average Length	242mm	(9.5")
	Xtra* 305mm	(12.0")

## Physical Properties

### Before Aging

Tensile Strength:	21 MPa
Ultimate Elongation:	550%

### After Aging

Tensile Strength:	21 MPa
Ultimate Elongation:	500%

## Quality Standards

Meets or exceeds ASTM D6319 - Standard Specification for Nitrile Examination Glove for Medical Applications.

Biocompatible per ISO 10993 studies for systemic toxicity, irritation and sensitization endpoints.

Passes Viral Penetration Testing - ASTM F1671 - 07 Standard Test Method for Resistance of Materials Used in Protective Clothing to Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage Penetration as a Test System.

## Clinical Considerations

Not formulated with Natural Rubber Latex – primary material is acrylonitrile butadiene, reducing the potential for glove-associated Type I allergic reactions to natural rubber latex proteins.

Powder-free, minimizing the potential for powder-related complications, such as irritant contact dermatitis.

**Tell your Halyard Health Sales Representative you want the durable feel and demonstrated performance of HALYARD\* PURPLE NITRILE\* Exam Gloves.**



**HALYARD  
Advantage\***

KNOWLEDGE NETWORK\* Accredited Education  
In-Service Training and Technical Support  
Credentialed Sales Representatives  
Tools & Best Practices  
Clinical Research  
Commitment to Excellence

For more information, please visit:

[www.halyardhealth.com](http://www.halyardhealth.com)

Call 1-844-HALYARD (1-844-425-9273)  
in the United States and Canada.

## HALYARD\* PURPLE NITRILE\* Exam Gloves

PURPLE NITRILE\* Exam Gloves have an AQL of 1.0 for critical defects. This means that PURPLE NITRILE\* Exam Gloves are significantly less likely to have critical defects than exam gloves that just meet the ASTM 2.5 AQL requirement.

### HALYARD\* PURPLE NITRILE\* Exam Gloves - 9.5"

Item	Description	Size	Packaging
55080	PURPLE NITRILE* Powder-Free Exam Gloves	X-Small	<b>100 eaches/box</b> 10 boxes/case
55081	PURPLE NITRILE* Powder-Free Exam Gloves	Small	<b>100 eaches/box</b> 10 boxes/case
55082	PURPLE NITRILE* Powder-Free Exam Gloves	Medium	<b>100 eaches/box</b> 10 boxes/case
55083	PURPLE NITRILE* Powder-Free Exam Gloves	Large	<b>100 eaches/box</b> 10 boxes/case
55084	PURPLE NITRILE* Powder-Free Exam Gloves	X-Large	<b>90 eaches/box</b> 10 boxes/case

### HALYARD\* PURPLE NITRILE-XTRA\* Exam Gloves - 12"

Item	Description	Size	Packaging
55090	PURPLE NITRILE-XTRA* Powder-Free Exam Gloves	X-Small	<b>50 eaches/box</b> 10 boxes/case
50601	PURPLE NITRILE-XTRA* Powder-Free Exam Gloves	Small	<b>50 eaches/box</b> 10 boxes/case
50602	PURPLE NITRILE-XTRA* Powder-Free Exam Gloves	Medium	<b>50 eaches/box</b> 10 boxes/case
50603	PURPLE NITRILE-XTRA* Powder-Free Exam Gloves	Large	<b>50 eaches/box</b> 10 boxes/case
50604	PURPLE NITRILE-XTRA* Powder-Free Exam Gloves	X-Large	<b>50 eaches/box</b> 10 boxes/case

### HALYARD\* PURPLE NITRILE\* Sterile Single Exam Gloves - 9.5"

Item	Description	Size	Packaging
52101	PURPLE NITRILE* Sterile Single Exam Gloves	Small	<b>100 eaches/box</b> 4 boxes/case
52102	PURPLE NITRILE* Sterile Single Exam Gloves	Medium	<b>100 eaches/box</b> 4 boxes/case
52103	PURPLE NITRILE* Sterile Single Exam Gloves	Large	<b>100 eaches/box</b> 4 boxes/case

### HALYARD\* PURPLE NITRILE-XTRA\* Sterile Pairs Exam Gloves - 12"

Item	Description	Size	Packaging
14260	PURPLE NITRILE-XTRA* Sterile Pairs Exam Gloves	Small	<b>50 eaches/box</b> 4 boxes/case
14261	PURPLE NITRILE-XTRA* Sterile Pairs Exam Gloves	Medium	<b>50 eaches/box</b> 4 boxes/case
14262	PURPLE NITRILE-XTRA* Sterile Pairs Exam Gloves	Large	<b>50 eaches/box</b> 4 boxes/case
14263	PURPLE NITRILE-XTRA* Sterile Pairs Exam Gloves	X-Large	<b>50 eaches/box</b> 4 boxes/case

[www.halyardhealth.com](http://www.halyardhealth.com)



\*Registered Trademark or Trademark of Halyard Health, Inc. or its affiliates. The COLOR PURPLE is a Registered Trademark of HYH or its affiliates. © 2015 HYH. All rights reserved.  
All other trademarks are the property of their respective owners.  
IP52 C14164

**Life-Assist**  
Helping Heroes Save Lives