

# BriteView® Color Code Laryngoscope Blades

Choose the BriteView® color code disposable blade system and know you have the right tool for the job

- **Super bright light** for maximum visibility compared to bent clear plastic rod design

- **Stainless steel from tip to heel** prevents flexing or breaking under stress



- **Single patient use** prevents cross contamination and eliminates cleaning costs
- **Color coded** for instant size recognition
- **Fits all standard handles** which are a fraction of the cost of fiber-optic types

## BriteView Color Code Blade System

Designed by EMS professionals with real-life street experience, the BriteView Laryngoscope Blade Series provides a strong, bright, laryngoscope blade, color coded for instant size recognition. Disposable stainless steel construction means you don't have to compromise strength and reliability for the cleaning and maintenance-free benefits of a single patient use blade.

The unique color coded heels make it a snap to select the correct size blade for the job at hand. The blade sizes are also etched on the side of the blade in a large format for easy selection when laid out on an instrument tray. The innovative package graphics are large and clear for easy identification in an emergency airway bag, under low light conditions or when restocking your vehicle's inventory.

If you want a disposable blade that won't flex, won't break off the handle and will provide a bright, reliable light source, then select the BriteView blade series today. Specify the BriteView color code blade on your next bid request. Contact your authorized Hartwell Medical dealer or give us a call at **1-800-633-5900**.

### HARTWELL MEDICAL

6354 Corte del Abeto, Suite F  
 Carlsbad, CA 92011  
 800-633-5900  
 760-438-5500 Tel  
 760-438-2783 Fax  
[www.HartwellMedical.com](http://www.HartwellMedical.com)

#### Available in Macintosh and Miller sizes



Model	Description (Sold by the box, 10 per box)
LB 42MIL0-DS	BriteView Miller 0 - White
LB 42MIL1-DS	BriteView Miller 1 - Pink
LB 42MIL2-DS	BriteView Miller 2 - Yellow
LB 42MIL3-DS	BriteView Miller 3 - Orange
LB 42MIL4-DS	BriteView Miller 4 - Lime Green
LB 42MAC1-DS	BriteView Macintosh 1 - Pink
LB 42MAC2-DS	BriteView Macintosh 2 - Yellow
LB 42MAC3-DS	BriteView Macintosh 3 - Orange
LB 42MAC4-DS	BriteView Macintosh 4 - Lime Green

*Creating Quality Products for Life*