## Non-Contact Infrared Body Thermometer





Measurement Time: ≤ 2 seconds

Measurement Distance: 0.4" – 2.0"

Body Temperature Range: 95.0 °F – 107.6 °F

Ambient Temperature: 59 °F – 95 °F

Accuracy: ± 0.4 °F

Automatic Power Off: ≤ 18 seconds

Batteries: 2 x AAA

Product Longevity: 3 years

Product Weight: 3 ounces (without batteries)

Product Size: L 5.43" x W 3.74" x H 1.47"

## This is an 510k cleared device

This non-contact infrared thermometer is intended to measure human body temperature on the forehead by consumers in a home setting and primary care settings as a screening tool.

## **MEDSOURCE**LABS

www.medsourcelabs.com

ltem #	Description	UOM
MS-131002	Non-Contact Infrared Thermometer, IR 300	40 EA/CS