



MedSource
Respiratory & Airway

Manual Suction Pump

MedSource Manual Suction Pump is a compact, hand powered, durable, light weight and high volume suction unit. Each unit is complete with an interchangeable canister, one adult and one child suction tube.



Features:

- Easy to operate by one hand
- Stable and shock resistant
- Heat resistant, reusable suction canister made of durable plastic that is shatter and leak proof
- Simple maintenance and easy to replace or clean canister
- Built-in overflow protection mechanism against spilling and clogging
- Reusable style pump
- Adjustable vacuum for child or adult
- Fits most common suction catheters
- User manual included

Specifications:

- Shipping Weight: 454 grams = 1 lbs
- Shipping Dimensions: 13x18x8cm = 9"x7"x3"

Rev. 1

www.gomedsource.com

Life-Assist 
Helping Heroes Save Lives

Item	Description	UOM
MS-001PMP	Manual Suction Pump	EA, 10 EA/CS
MS-002PMP	Refill Canister	EA, 50 EA/CS