

Protective Gown #86792

Disposable Poly Protective Gown

- Long sleeve, thumb loop.
- Over-the-head, easy-tear neck.
- Waist tie and open back.
- AAMI Level 3.
- Unisex.
- One size fits most.
- Greatly reduce the risk of cross-contamination.
- Enjoy the comfort & confidence from a fresh, new and clean product every time.



Item #	Description	Material	Color	Size	Case	Quality	AAMI Level
86792	Protective Gown	Poly	Blue	42" x 46"	1/BG; 75BG/CS	Economy	Level 3