

NEW

REX

RapidEXtraction

Rescue at the Speed of Life when seconds count



The Rex Rapid Extraction Mobile Stretcher's rugged and durable design provides rapid patient movement in the roughest conditions. The REX's multi-functional design does the work of a stretcher, gurney, backboard and decent chair all with one system.

REX is the only patient transport product that provides resource advantages allowing a single operator to move injured faster than traditional transport devices.



SWAT Active Shooter



Mass Casualty Treatment Cot



Fire Equipment Support



EMS Events & Venues

When Time is Critical, Lives Depend on *Rapid Extraction*



Patent Number
US 7,549,661 B1

REX Rapid Extraction Platform designed for rapid deployment by a single responder. The 360° maneuverability allows responders to safely move patients of any size.

REXONE



Designed for emergency extractions and evacuations in any environment.

- Rugged T1 Aerospace aluminum frame
- 42 lbs. fully assembled
- 75"x25" fully assembled-folded: 26"x18"
- Height from ground to base 10"
- 400 lbs. weight capacity
- Tri-fold panel system for quick assembly and storage. Converts from backboard, mobile stretcher, to patient treatment cot within seconds
- Quick release wheel system for interchangeable wheels
- Telescoping pull handle
- Adjustable locking head rest
- Interchangeable patient straps (set of 3)
- Backpack carrying case

REX Descent Control System



Provides all RexOne functionality, specifications, and accessories. Evacuates down stairwells or steep terrains.

- Duel rod braking system for even braking in any condition
- Hand Brake Control for ease of use. Controls speed and provides braking ability during descent
- Emergency Brake locking System prevents rolling during decent downstairs or steep terrain
- 2 point Safety Connection Brackets and safety straps provides additional support for responders during difficult descents
- Adjustable locking head rest
- Interchangeable patient straps (set of 3)
- Backpack carrying case

Accessories:

ALL TERRAIN TIRES	SAND/ICE/SNOW TIRES	STRAP SYSTEM	ADJUSTABLE IV POLE	BIKE TOW SYSTEM	BACKPACK CARRYING CASE
 <p>20" puncture free for all conditions. Handles rough conditions and requires no maintenance. Weather resistant from heat and cold.</p>	 <p>Sand- Designed for beach and sandy condition.</p>  <p>Snow/Ice Tires- Handles light snow and icy condition as well.</p>	 <p>Black – set of 3 Tri-Colored – set of 3 special order Patient Strap Extenders—Provides additional length for chest and mid-section straps</p>	 <p>Removable pole for easy storage</p>  <p>Safety Bracket Safety Support Straps - provides additional for support responders. Used to secure REX during vehicle transport</p>	 <p>Attached to most street police bike frames</p>	 <p>Adjustable straps with external wheel travel pocket for added protection</p>

