

# Technical Data Sheet

## **Polypropylene Metal Double D-Ring 20' Strap**



Item # 250201

**Contact Information:**

Distributor: Dick Medical Supply, LLC.  
Address: 1946 Camaro Ave,  
Columbus, Ohio 43207, USA  
Phone: (614)-444-2300  
Fax: (614)-444-2301  
E-Mail: [kfannin@dms4ems.com](mailto:kfannin@dms4ems.com)  
Revision: 02152022



## 1. Product Overview

- Designed to minimize patient movement
- Incorporates 2” Polypropylene Webbing, 2” Metal Double D-Ring
- Strap is fully adjustable for a custom fit
- Fits standard Stryker and Ferno Ambulance Cots

## 2. Specifications

Product Name:	Polypropylene, Metal Double D-Ring Buckle, 20 ft., 1 piece strap
Color:	Available in Black (BK), Blue (BL), Green (GR), Neon Green (NG), Orange (OR), Red (RD), White (WH), Yellow (YL)
Size:	2” Wide X 20’ Long
Certification:	CE, FDA
Reusable:	Yes
Country of Origin:	USA
Application:	EMS
Item#:	250201
Pack Size:	1ea

**Disclaimer:** Dick Medical Supply, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. Dick Medical Supply makes no representations or warranties either express or implied, regarding the suitability of the material for any purpose or the accuracy of the information contained within this document. Accordingly, Dick Medical Supply will not be responsible for damages resulting from use of or reliance upon this information.