



RÜSCH ASHERMAN CHEST SEAL™

Quick, efficient and safe treatment of open chest injuries

AT A GLANCE

Open chest wound dressing set for effective emergency management:

- unique one-way valve
- transparent circular pad, 14 cm in diameter
- strong adhesive
- includes a gauze pad
- sterile
- easy-to-open package

The Teleflex RÜSCH brand Asherman Chest Seal™ (ACST™) is a safe approach to treating open pneumothorax and preventing tension pneumothorax in chest injuries resulting from gunshot, stab wounds, and other penetrating chest trauma.

Its unique one-way valve allows air & blood to escape, while preventing their re-entry.

The sterile chest wound dressing set has a number of distinguishing features, making it highly effective in cases of emergency:

- a unique one-way valve for maximised safety
- an ample dressing size of 14 cm in diameter
- a circular pad, made of a transparent material that allows optimal visualisation of the wound
- a gauze pad for cleaning and drying the wound area

- a pressure-sensitive adhesive that ensures an effective seal, even over body hair

Proven to be more effective than standard petroleum gauze, the ACST™ is standard issue for the U.S. Army, the British Army and the U.S. Navy.

Ordering information

ASHERMAN CHEST SEAL RÜSCH	
REF.	QTY
849100	10

DISTRIBUTED BY:

TELEFLEX MEDICAL HEADQUARTER EMEA, IRELAND · IDA Business Park · Athlone · Co Westmeath

Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com

UNITED KINGDOM Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com

WWW.TELEFLEX.COM

94 03 66-00 00 01 · REV A · 08 12 01 · All data current at time of printing. Subject to technical changes without further notice.