



**MedSource**  
Infection Control

# Sharps Dart



**510K Cleared**

The MedSource Sharps Dart features a simple, one-handed operation to open and close the disposal opening. The convenient size fits easily into carry bags for EMS, Fire/Rescue and Nursing for containment of used IV catheters and smaller syringes.

The leak-resistant cap offers a secure closure between uses when properly closed and the new locking lid ensures containment of contents. The large disposal opening easily accommodates small syringes and other sharps waste.

Rev. 1

[www.gomedsource.com](http://www.gomedsource.com)

Item #	Description	UOM
MS-64250	Sharps Dart	24 EA/CS, 8 CS/MC, 192 EA/MC

## Features:

- Small Size
- One handed Operation
- Locking Lid
- Puncture resistant bottom

## Benefits:

- Portable for EMS Professionals
- Safely dispose of sharps immediately after use.
- Once sharp is contained with a locked lid risk of accidental needlestick is eliminated.
- The heavy weight plastic bottom has been tested and proved puncture resistant.

## Specification: MS-64250

### Materials

Plastic

### Dimensions

Length: 6.7" (17cm)  
Diameter: 1.2" (3cm)  
Empty Weight: .9oz (255grams)

### Product Testing

Tested to the Standard of ISO ASTM F2132-01 for Puncture Resistance