

TECHNICAL DATA SHEET


Elastic Bandage with Velcro

Intended Use

Elastic bandage with Velcro used to provide flexible support and compression during activity. Ideal for strains and sprains. This type of product fits all parts of the body. To be used in health care facilities or over the counter single use.

Characteristic	Acceptable Limit		
Material Composition	6 - 16% Spandex, 35 - 45% Polyester, 44 – 54% Polyester-Cotton. Latex-free.		
Description	Roll of woven elastic material with Velcro on the core end. The end of Velcro and the coil are stitched with a color (Refer to Appendix 1 specific attribute)		
Dimension	Size	Width (cm)	Stretched Length (m)
	2" x 5 YD	5,1	≥ 4,2
	3" x 5 YD	7,6	
	4" x 5 YD	10,2	
	6" x 5 YD	15,2	

Visual



Additional Information

Storage Condition	Country of Origin
Store under ambient temperature and dry conditions	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment
A632	AMD-Ritmed Elastic Bandage Roll with Velcro	1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case Individually Wrapped Wrap has a perforated line circling the roll in the middle of the length	2" x 5 YD (5,0 cm x 4,5 m)	White and skin color stripes	The Velcro and the end of the coil are stitched with yellow color	Non-Sterile
A633	AMD-Ritmed Elastic Bandage Roll with Velcro	1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case Individually Wrapped Wrap has a perforated line circling the roll in the middle of the length	3" x 5YD (7,5 cm x 4,5 m)	White and skin color stripes	The Velcro and the end of the coil are stitched with red color	Non-Sterile
A634	AMD-Ritmed Elastic Bandage Roll with Velcro	1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case Individually Wrapped Wrap has a perforated line circling the roll in the middle of the length	4" x 5YD (10,0 cm x 4,5 m)	White and skin color stripes	The Velcro and the end of the coil are stitched with blue color	Non-Sterile
A636	AMD-Ritmed Elastic Bandage Roll with Velcro	1 Roll/Wrap, 10 Wraps/Box, 5 Boxes/Case, 50 Rolls/Case Individually Wrapped Wrap has a perforated line circling the roll in the middle of the length	6" x 5 YD (15,0 cm x 4,5 m)	White and skin color stripes	The Velcro and the end of the coil are stitched with green color	Non-Sterile