



Burn Kits

Be Prepared for Any Burn Emergency.

More than 30 years ago, Water-Jel Technologies introduced water-based gel technology that revolutionized emergency first response for burn injuries. Today, WATER-JEL is used by Emergency, Fire-Rescue, Industrial Safety and Military professionals around the world.

When Seconds Count

The “immediate cooling” of burns is recommended by EMS protocols and curriculum. Water-Jel Burn Dressings cool burns on contact to relieve pain fast while covering and protecting burn injuries. They consist of a medical-grade non-woven pad, pre-moistened with Water-Jel’s NEW clear, fragrance-free, water-based gel with Hyaluronic Acid. The gel is made from 96% water and is slow to evaporate, so blankets and dressings remain moist and non-adherent. They provide thermal cooling of burns (not evaporative) without the need for re-wetting so they won’t get too cold or increase the risk of hypothermia – even on large extent burns. Safe and effective, Water-jel burn care is easy to use for burn emergencies, when seconds count.

First Response, Safety and Preparedness

Water-Jel Burn Kits come in different sizes, featuring unique product assortments for different work environments.

- Fire rescue
- EMS
- Manufacturing
- Electric/gas utilities
- Petrochemical
- Foundries
- Law enforcement
- Military

Specifications

- Pre-moistened, easy to use, won’t stick to burn injuries.
- Available in 6 sizes: 2” x 6”, 4” x 4” , 4” x 16”, 8” x 18”, 8”x20”, 12”x16”.
- Special burn dressings for Hand Burns and Facial Burns.
- Medical Device - contains no drugs/active ingredients.
- Made in the USA at an FDA-regulated facility.

Soft-sided burn kits can be stored on a wall and also in vehicles behind and under a seat. Hard plastic Burn Kits and Burn Stations come with mounting holes.



WATER-JEL BURN KITS



LARGE SOFT-SIDED BURN KIT WITH BURN WRAP ITEM EBK2-HA

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 One 8" x 18" Burn Dressing
 One 4" x 16" Burn Dressing
 Three 4" x 4" Burn Dressing
 Four sterile gauze roll 3" x 4 yds.
 One 36" x 30" (3' x 2.5') Fire Blanket/Burn Wrap
 One scissors
 One instructions card
 Heavy-duty nylon carry bag with inner pockets.

- Vital in areas where high heat flames, arc flash, hot liquid, steam or chemicals exist.
- Everything needed for Industrial, Institutional and Medical burn emergencies.

Case Pack: 3 Kits
Case Size: 15.38" x 12.63" x 8.13"
Case Weight: 22.3 lbs
Case Cube: 0.92 cu. ft.



SMALL SOFT-SIDED BURN KIT ITEM EBK1-HA

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 Two 4" x 4" Burn Dressing
 Two 4" x 16" Burn Dressing
 Two sterile gauze roll 3" x 4 yds.
 One scissors
 One instructions card
 Heavy-duty nylon carry bag with inner pockets.

- Perfect for use in work vehicles.
- Durable nylon handle straps for easy carrying, or can be hung on a wall.

Case Pack: 6 Kits
Case Size: 15.38" x 12.63" x 8.13"
Case Weight: 13.7 lbs
Case Cube: 0.92 cu. ft.



SOFT-SIDED CRITICAL BURN KIT ITEM BKCK-HA

Two 4" x 4" Burn Dressing
 One 8" x 20" Hand Burn Dressing
 One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 Two sterile gauze rolls 3" x 4 yds., One pair of scissors

- Carry Case: High-impact polystyrene kit box has living hinges, strong nylon catches, moisture-resistant cover gasket and fold-down combination handle/hanger.

Case Pack: 6 Kits
Case Size: 15.38" x 12.63" x 8.13"
Case Weight: 19.9 lbs
Case Cube: 0.92 cu. ft.



INDUSTRIAL/WELDING BURN KIT ITEM IWK-HA

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 One 8" x 18" Burn Dressing
 One 4" x 16" Burn Dressing
 Two 4" x 4" Burn Dressing
 Two 2" x 6" Burn Dressing
 One BurnJel dispenser with 25 unit dose packets
 One instructions card • One elastic bandage

- Carry Case: High-impact polystyrene kit box has living hinges, strong nylon catches, moisture-resistant cover gasket and fold-down combination handle/hanger. Can be wall-mounted in locations where burns can happen.

Case Pack: 5 Kits
Case Size: 16" x 12" x 8"
Case Weight: 23.3 lbs
Case Cube: 0.89 cu. ft.



FOOD SERVICE BURN KIT ITEM FSK-HA

One 12" x 16" Burn Dressing (with pre-cut holes for facial burns)
 Two 4" x 4" Burn Dressing
 Two 2" x 6" Burn Dressing
 One BurnJel dispenser with 25 unit dose packet
 One instructions card • One elastic bandage

- Carry Case: High-impact polystyrene kit box has living hinges, strong nylon catches, moisture-resistant cover gasket and fold-down combination handle/hanger. Can be wall-mounted in locations where burns can happen.

Case Pack: 5 Kits
Case Size: 16" x 12" x 8"
Case Weight: 12.3 lbs
Case Cube: 0.89 cu. ft.



LARGE EMERGENCY BURN STATION ITEM EBSL-HA

One Face Dressing 12" x 16"
 One Dressing 8" x 18"
 Two Dressings 4" x 16"
 Two Dressings 4" x 4"
 One Burn Jel 2 oz bottle
 Ten Burn Jel unit dose 3.5g
 One pair of scissors
 Two rolls of gauze

- Provides immediate access to a wide range of burn treatment products.
- Place in strategic locations where burn accidents are most likely to occur.
- Comes with mounting holes and hardware.

Case Pack: 5 stations
Case Size: 17.63" x 22" x 12.75"
Case Weight: 12.55 lbs
Single Pack Size:
 17.25" x 21.25" x 2.25"
Weight: 2.33 lbs



SMALL EMERGENCY BURN STATION ITEM EBSS-HA

Two Dressings 2" x 6"
 Two Dressings 4" x 4"
 Six Burn Jel unit dose 3.5g
 One roll of gauze

- Offers full product visibility.
- Very well suited for fast food operations and smaller businesses.
- Comes with mounting holes and hardware.

Case Pack: 5 stations
Case Size: 15.5" x 8.75" x 9.25"
Case Weight: 2.73 lbs
Single Pack Size:
 14.875" x 8.5" x 1.5"
Weight: .47 lbs



COMPACT BURN KIT ITEM BK10-HA

One Dressing 2" x 6"
 One Dressing 4" x 4"
 One sterile gauze roll 3" x 4 yds.
 Six Burn Jel unit dose 3.5g
 One pair of scissors

Case Pack: 5 Kits
Case Size: 16" x 12" x 8"
Case Weight: 23.3 lbs
Case Cube: 0.89 cu. ft.



COMPACT BURN KIT PLUS ITEM BK11-HA

One Dressing 2" x 6"
 One Dressing 4" x 4"
 One Dressing 4" x 16"
 Six Burn Jel unit dose 3.5g
 One pair of scissors
 One sterile gauze roll 3" x 4 yds.

Case Pack: 5 Kits
Case Size: 16" x 12" x 8"
Case Weight: 23.3 lbs
Case Cube: 0.89 cu. ft.