

SANDEL® Z-SLIDER™ PATIENT TRANSFER SHEET


STAFF & PATIENT ERGONOMICS

The #1 Best-Selling Disposable Patient Transfer Sheet

The Z-Slider™ disposable patient transfer sheet helps prevent disabling back injuries that can result from lifting, transferring and repositioning patients. Z-Slider™ is placed under the draw sheet and acts to reduce friction while moving the patient.


Use of Z-Slider™ complies with recommendations by OSHA and the AORN for safe patient handling, as well as “No-Lift” laws in effect in many states.^{1,2}

To enhance staff compliance, Ansell offers a variety of point-of-care dispensers.



KEY FEATURES & BENEFITS

- NEW! Now available in XL size to accommodate tall and bariatric patients
- Reduces the risk of physical strain to the back, shoulders, neck and arms caused by moving patients
- Disposable to help prevent cross-contamination
- No weight limit
- Cost-effective
- Easy to use - product instructions are printed on the product
- Recyclable (LDPE 4)
- Standard Size: 33” x 39”
XL Size: 45” x 53”



SANDEL® Anticipate the Unexpected

Developed in collaboration with clinical professionals, SANDEL products are designed to address key high-risk areas in the healthcare environment. Our signature bright orange products help reduce the risk of preventable injuries, errors and infections. SANDEL products equip caregivers with solutions to do their job safely while improving patient outcomes.

¹ Hospital eTool: HealthCare Wide Hazards Module: Ergonomics, U.S. Department of Labor, Occupational Safety and Health Administration (OSHA), www.osha.gov
² Pocket Reference Guide: Safe Patient Handling and Movement in the Perioperative Setting, AORN Toolkit, Association of periOperative Nurses, 2014

SANDEL® Z-SLIDER™ PATIENT TRANSFER SHEET

STAFF & PATIENT ERGONOMICS

ITEM #	PRODUCT NAME	DESCRIPTION	PACKAGING	CASE QTY.
2101	Z-Slider™ Single Pack	Patient Transfer Sheet, Individually Wrapped, 33" x 39" (Non-Sterile)	24/Bx, 2 Bx/Cs	48
2151	Z-Slider™ XL Single Pack	Patient Transfer Sheet, Individually Wrapped, 45" x 53" (Non-Sterile)	12/Bx, 2 Bx/Cs	24
2101-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Dispenser for use with Z-Slider™ and Z-Slider™ XL Single Packs	4/Cs	4
2102	Z-Slider™ - Roll	33" x 39" Patient Transfer Sheet in Roll Dispenser (Non-Sterile)	30/Bx, 2 Bx/Cs	60
2103	Z-Slider™ Wall Dispenser (Large)	Patient Transfer Sheet Dispenser Refill for Large Point-of-Care Wall Holder (33" x 39" Patient Transfer Sheet, Non-Sterile)	15/Bx, 4 Bx/Cs	60
2103-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Mounting Holder for the 2103 Dispenser Box (12" x 20.5")	5/Cs	5
2104	Z-Slider™ Wall Dispenser (Small)	Patient Transfer Sheet Dispenser Refill for Small Point-of-Care Clear Wall Holder (33" x 39" Patient Transfer Sheet, Non-Sterile)	10/Bx, 6 Bx/Cs	60
2104-H	Point-of-Care Clear Wall Holder	Point-of-Care Clear Wall Mounting Holder for the 2104 Dispenser Box (11.5" x 11")	5/Cs	5



2101
Z-Slider™ Single Pack
Individually Wrapped



2151
Z-Slider™ XL Single Pack
Individually Wrapped



2101-H
Z-Slider™ Single Pack
Point-of-Care Clear
Wall Mounting Holder



2102
Z-Slider™ Roll



2103
Z-Slider™ Refill for Large
Point-of-Care Clear
Mounting Wall Holder



2103-H
Z-Slider™ Large Point-of-Care
Clear Wall Mounting Holder



2104
Z-Slider™ Refill for Small
Point-of-Care Clear Wall Holder



2104-H
Z-Slider™ Small Point-of-Care
Clear Wall Mounting Holder

➤ **Contact your Ansell representative for ordering or more information.**

North America
US Tel: 866-764-3327
818-534-2500
Email: sandel-cs@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0)2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com