



# PNEUMERIC CAPNOSPOT

Decompress the Chest with Objectivity and Confidence!

## Capnospot – First Pass Success

- Decompression of tension pneumothorax via needle thoracostomy is lifesaving
- Prehospital Decompression Reduces 24 hour Mortality by 25%<sup>(1)</sup>
- The current method of judging needle decompression success using a “gush of air” and clinical signs is subjective.
- **Even well trained providers with correct equipment have Needle Thoracostomy failure rates of 20-50%**<sup>(2,3)</sup>



- Capnospot provides objective visual confirmation of decompression success or failure by utilizing visual color change, ideal for noisy prehospital environments.
- Small, Lightweight, Compact
- Capnospot is compatible with all decompression devices.
- The Capnospot is more accurate than the standard of care<sup>(4)</sup>
- Potential to improve decompression rates, reduce the number of procedures, reduce wasted procedures and costs. <sup>(4)</sup>
- Objective color change is nearly instant (<5 seconds) <sup>(4)</sup>
- Capnospot indicates if the decompression is non-therapeutic, or catheter loses patency will revert to its original color. <sup>(4)</sup>
- Functions in Low light environments



## Capnospot – **KNOW** You are in!



[www.pneumeric-medical.com](http://www.pneumeric-medical.com)

**Life-Assist**  
Helping Heroes Save Lives

1. Muchnok D, Vargo A, Deeb AP, Guyette FX, Brown JB. Association of Prehospital Needle Decompression With Mortality Among Injured Patients Requiring Emergency Chest Decompression. JAMA Surg. 2022 Oct 1;157(10):934-940. doi: 10.1001/jamasurg.2022.3552. PMID: 35976642; PMCID: PMC9386601.  
2. J. M. Aho, C. A. Thiels, M. M. El Khatib, D. S. Ubl, D. V. Laan, K. S. Berns, E. B. Habermann, S. P. Zietlow, M. D. Zielinski, Needle thoracostomy: Clinical effectiveness is improved using a longer angiocatheter. J. Trauma Acute Care Surg. 80, 272–7 (2016).  
3. D. V. Laan, T. D. N. Vu, C. A. Thiels, T. K. Pandian, H. J. Schiller, M. H. Murad, J. M. Aho, Chest wall thickness and decompression failure: A systematic review and meta-analysis comparing anatomic locations in needle thoracostomy. Injury 47, 797–804 (2016).  
4. Naik ND, Hernandez MC, Anderson JR, Ross EK, Zielinski MD, Aho JM. Needle decompression of tension pneumothorax with colorimetric capnography. Chest. 2017;152:1015–20. - DOI