

CaviWipes™ HP

Surface Disinfectant

CaviWipes HP disinfecting wipes are the newest members of the trusted CaviWipes family providing a 1-minute universal contact time from a hydrogen peroxide formulation that qualifies for the EPA's rigorous Emerging Viral Pathogen Claim². The High Performance disinfection protocol offered by CaviWipes HP helps protect both patients and staff.



Efficient Disinfection

1-Step cleaning/disinfection³ and 1-minute universal contact time for all listed pathogens ensures efficient disinfection.



Broad Spectrum Efficacy

In addition to providing effective disinfection against all 67 listed pathogens including MDRO's (Multidrug Resistant Organisms), CaviWipes HP disinfecting wipes qualify for the EPA's rigorous Emerging Viral Pathogen claim², protecting you against future pathogens as well.



Ease of Use

The new green top opens and closes more easily and securely, and dispenses wipes with precision, helping to reduce waste.



Lowest Toxicity Level

US EPA Toxicity Category IV products minimize exposure to harmful cleaning chemicals. The use of gloves is not required.⁵



1 Minute Efficacy Against:⁴

Mycobacterium

Mycobacterium tuberculosis var: *bovis* (BCG)

Bacteria

Pseudomonas aeruginosa
Staphylococcus aureus
Salmonella enterica
Bordetella pertussis
Enterococcus faecium
Klebsiella aerogenes
Klebsiella pneumoniae
Escherichia coli
Escherichia coli O157:H7
Burkholderia cepacia
Stenotrophomonas maltophilia
Streptococcus pyogenes
Listeria monocytogenes
Shigella dysenteriae
Serratia marcescens
Bordetella bronchiseptica
Campylobacter jejuni
Proteus vulgaris

Multi-Drug-Resistant Bacteria (MDRO)

Penicillin-Resistant *Streptococcus pneumoniae* (PRSP)
 Vancomycin-Resistant *Enterococcus faecalis* (VRE)
 Methicillin-Resistant *Staphylococcus aureus* (MRSA)
 Hospital-Acquired and Community Acquired Methicillin-Resistant *Staphylococcus aureus* (HA-MRSA and CA-MRSA)
 Vancomycin Intermediate *Staphylococcus aureus* (VISA)
 Vancomycin Resistant *Staphylococcus aureus* (VRSA)
 Linezolid-Resistant *Staphylococcus aureus* (LRSA)
 Carbapenem-Resistant *Klebsiella pneumoniae*
 Multidrug-Resistant *Acinetobacter baumannii* (MRAB) Strains ATCC 19606 and ATCC BAA-1605
 NDM-1 *Enterobacter cloacae*
 NDM-1 *Klebsiella pneumoniae*
 MBL *Pseudomonas aeruginosa*
 NMD-1 *Escherichia coli*
 ESBL *Klebsiella pneumoniae*
 ESBL *Escherichia coli*
 Multidrug-Resistant *Klebsiella pneumoniae*
 Carbapenem-Resistant *Escherichia coli*

Pathogenic Fungi

Trichophyton interdigitale
Candida auris
Candida albicans

Viruses (Non-Enveloped)

Human Norovirus
 Human Rotavirus
 Rhinovirus Type 37
 Rhinovirus Type 1A
 Human Adenovirus Type 5
 Echovirus Type 11
 Enterovirus Type 71

Viruses (Enveloped)

Herpes simplex virus Type 1 & Type 2
 Avian Influenza A virus (H7N9)
 Avian Influenza A virus (H5N1)
 Influenza A virus (H3N2)
 Influenza B virus (H1N1)
 Parainfluenza virus Type 3
 Hepatitis B virus (HBV)
 Human Hepatitis C virus (HCV)
 Human Immunodeficiency virus Type 1 (HIV-1)
 Severe Acute Respiratory Syndrome-related Coronavirus 2 (SARS-CoV-2)
 Vaccinia virus
 SARS-Associated Coronavirus (SARs-CoV)
 Cytomegalovirus (CMV)
 Influenza A virus (H1N1) Strains A/4108/2009 and A/04/09
 Human Coronavirus
 Respiratory Syncytial Virus (RSV)
 Measles virus

Active Ingredient

Hydrogen Peroxide	0.51%
Other Ingredients†	99.49%
Total	100.00%

[† Does not include the weight of the towelette]



Ordering Information

Part No. Description

16-1100	CaviWipes HP (6" x 6.75") – 160 towelettes per canister
16-1150	CaviWipes HP XL (9" x 12") – 65 towelettes per canister

- List N: Disinfectants for Use Against SARS-CoV-2. Feb 28, 2022, www.epa.gov/pesticide-registration/list-n-disinfectants-use-against-sars-cov-2.
- CaviWipes HP qualifies for the EPA Emerging Viral Pathogen program against all virus tiers (enveloped, small and large non-enveloped).
- Follow Product Label for special precleaning and PPE instructions.
- Always follow Product Label for Cleaning/ Disinfection Directions.

To learn more about CaviWipes HP scan the QR code or visit metrex.com

