



Defender™ DFR10

Nitrile Exam Gloves



Nitrile

Meets the unique demands of First Responders

- Advanced non-latex nitrile polymer offers superior stretch and comfort
- Textured fingertips for enhanced grip
- AQL of 0.65 exceeding industry standards
- Reduced hand fatigue & improved dexterity
- No-stick exterior surface prevents tape from sticking to the glove




Industries

- Emergency Medical Services
- Healthcare
- Law Enforcement

Recommended For

- Administering drugs
- Standard patient examinations
- Emergency services





Defender™ DFR10

Nitrile Exam Gloves

TECHNICAL DATA SHEET

PRODUCT INFORMATION

	Defender™ DFR10
Material	Nitrile
Color	Violet Blue
Glove Design	Powder-Free, Chlorinated, Textured Fingers
Cuff	Beaded
Manufacturing/QMS Audit Standards	ISO 13485, ISO 9001
Regulatory/Standards Compliance	ASTM D6319, FDA21 CFR 177-2600, ISO 13485
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	Malaysia
User Needs Segment	Robust
Available sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10), XXL (10.5 - 11)
Anti-static	Not Tested

PHYSICAL PROPERTIES

	Typical Values		Testing Method
Length (mm/inches)	250 / 9.8		ASTM D6319
Freedom from Holes (Inspection level I)	0.65 AQL		ASTM D5151
Palm Thickness (mm/mils)	0.13 / 5.1		ASTM D6319
Finger Thickness (mm/mils)	0.19 / 7.5		ASTM D6319
	BEFORE AGING	AFTER AGING	
Ultimate Tensile Strength (MPa)	24	20	ASTM D412 & D573
Elongation at Break (%)	670	505	ASTM D412 & D573

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)	XXL (10.5 - 11)
Product Code	DFR10-1	DFR10-2	DFR10-3	DFR10-4	DFR10-5

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa Region

Ansell Healthcare Europe NV
T: +32 (0) 2 528 74 00
F: +32 (0) 2 528 74 01

Asia Pacific Region

Ansell Global Trading Center
(Malaysia) Sdn Bhd
T: +603 8310 6688
F: +603 8310 6699

North America Region

Ansell Healthcare Products LLC
T: +1 800 800 0444
F: +1 800 800 0445

Latin America & Caribbean Region

Ansell Commercial Mexico S.A. de C.V.
T: +52 442 248 1544 / 248 3133

Australia

Ansell Limited
T: +61 1800 337 041
F: +61 1800 803 578

Russia

Ansell PYC
Ten. +7 495 258 13 16

Performance Standards and Regulatory Compliance



ISO 13485

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: www.ansell.com/patentmarking © 2021 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

