

# FASPLINT®

## Extremity Splint Kit



### Item Number: SP912C

#### All-Aluminum COMPACT Vacuum Pump (3 Year Warranty)

**Pump dimensions:**

- 11.0" Length (Top of handle to bottom of pump cylinder)
- 11.0" Length for storage with hose connected
- 2.5" Width (pump cylinder only)
- 2.5" Width (pump cylinder and end caps)
- 4.0" Handle width

**Pump weight:** 2.2 lbs.

**Pump materials:** Aluminum, stainless, brass and acetal

**Pump Shaft:** Stainless

**Interior O-ring:** Rubber, replaceable

**Hose Length:** 42"

**Hose Length (Cylinder to end of hose end connector):** 42.8"

**Hose Material:** Dual wall, vinyl exterior, rubber interior

**Hose end connector:** Plastic tapered tip

**Hose connection point:** Recessed in side of pump cylinder



### Item Number: SP912

#### Economy Vacuum Pump (6 Month Warranty)

**Pump dimensions:**

- 12.0" Length (Top of handle to bottom of pump cylinder)
- 13.5" Length for storage with hose connected
- 2.4" Width (pump cylinder only)
- 3.0" Width (pump cylinder and end caps)
- 4.0" Handle width

**Pump weight:** 0.8 lbs.

**Pump materials:** Plastic and aluminum

**Pump Shaft:** Plastic

**Interior O-ring:** None, Plastic Flange

**Hose Length:** 27"

**Hose Length (Cylinder to end of hose end connector):** 28.5"

**Hose Material:** Single Wall Silicone

**Hose end connector:** Plastic tapered tip

**Hose connection point:** Exposed at bottom of Pump

DISTRIBUTED BY 

**Life-Assist**  
Helping Heroes Save Lives 