


SIMBODIES®

Realistic training. Real **life-saving.**

EMERGENCY MEDICAL SKILL MANIKIN (EMS-MAN)



Improve your lifesaving skills on the world's most hyper-realistic adult medical manikin. Thoughtfully designed and sculpted from real human anatomy by clinical experts, scientists, and artists.

Key Features

Life-like look and feel with flexible limbs

Partially open eyes

Accurate human cast cadaver airway

Age: Approx. 30-40 years old

Skin tone - White Caucasian

Clothing size: Medium/Large

Shoe size: UK 10, EU 45, US 11

Height: 6' Weight: 35kg/ 77lbs

Includes:

- 2 x Amputation Sleeves
- 2 x Manual Blood Pumps

Maintenance kit:

- 100ml Bottle of Sweat & Tears (lubrication)
- Powder and Brush
- Bluetooth speaker, manual, and charging cable

** additional items included*

SIMBODIES®
Realistic training. Real **life-saving.**
Simbodies.com
SafeguardMedical.com

DISTRIBUTED BY 
Life-Assist
Helping Heroes Save Lives 

To order from the SIMBODIES line of products today or for more information, talk with your representative or contact Safeguard Medical at 855-428-6074. www.safeguardmedical.com

PART OF **SAFEGUARDMEDICAL™**

SIMBODIES®

Realistic training. Real **life-saving.**

MEDICAL PROCEDURE	ADULT MALE
Catastrophic/ Life-Threatening Hemorrhage	✓
Application of a Tourniquet	✓
Wound Packing	✓
Junctional Bleeding	✓
Airway Management	✓
Head Tilt	✓
Chin Lift	✓
Jaw Thrust	✓
Postural Airway Management (1)	✓
Nasopharyngeal Airway Insertion (1)	✓
Oropharyngeal Airway Insertion (1)	✓
Supraglottic Airway Device Insertion (2)	✓
Laryngoscopy/ Tracheal Intubation (1)	✓
Needle Cricothyroidotomy (9)	✓

1. For most current medical devices 2. Different supraglottic airway devices will have varying degrees of realism in insertion and ventilation, depending on the device used. 3. Chest expansion during ventilation can vary depending on the manikin, airway, and ventilation devices used 4. Some chest seals may not adhere to silicone skin in the same manner as human skin. 5. Child/baby and Adult Female 2nd intercostal only, Adult Male 2nd and 5th 6. Chest compressions can be performed on all of our manikins but are not designed solely for CPR training. Chest Compressions depth and feedback varies by manikin. 7. Certain devices will work. Please contact simbodiesupport-us@safeguardmedical.com 8. Some pads may not adhere to silicone skin in the same manner as human skin. 9. Available on the EMS-T Trainer.

C-Spine Control	✓
Application of a Spinal Collar	✓
Application of Manual Inline Stabilisation	✓
Breathing Assessment	✓
Artificial Ventilation (BVM) (3)	✓
Artificial Ventilation (Ventilator when intubated) (3)	✓
Application of Chest Seal (4)	✓
Needle Thoracocentesis (5)	✓
Circulation Assessment	✓
Application of Various Dressings	✓
Application of Pelvic Binder	✓
Application of Traction Device (Femoral)	✓
Intraosseous Access Tibia	✓
Intraosseous Access Humeral Head	✓
Chest Compressions (manual) (6)	✓
Chest compressions (automated) (7)	✓
Application of AED/ Defibrillator Pads (8)	✓

* no customizable options available on this product



SIMBODIES®
Realistic training. Real **life-saving.**

Simbodies.com
SafeguardMedical.com

DISTRIBUTED BY 
Life-Assist
Helping Heroes Save Lives 

To order from the SIMBODIES line of products today or for more information, talk with your regional representative or contact simbodies.com/contact-us

PART OF **SAFEGUARDMEDICAL**