



The G3 RESPONDER is a large-sized street pack for medics who carry more. Boxy and compact with ergonomic hauling.

Unique design keeps your shoulder straps off the ground when you're accessing contents. Internal Cells stack neatly and securely in many different configurations. This pack accommodates the First Aid QuickRoll Intubation Kit, Circulatory Kit, and Remedy Kit, as well as the stackable Cellular System. The aluminum track creates a bomb proof secure point for the shoulder straps that allows for comfortable transport of most emergency medical equipment. Up to 4 of the Generation 3 standard Cells will fit inside as this is the Largest size of the EMS Backpacks that are designed to use the Cellular System.

Dimensions: 28.5" x 17" x 7"

Weight: 6 lbs

Available colors: ●

SKU: G35000RE



Learn more: statpacks.com/g3-responder