

# VIROMAX™ IN-LINE BVF

BIDIRECTIONAL BACTERIAL-VIRAL FILTER FOR IN-LINE USE

**ViroMax™ bacterial-viral filter features coaxial ISO connectors for applications in many stages of breathing circuits and mechanical ventilation.**

- Bacterial filtration efficiency: > 99.999%
- Viral filtration efficiency: > 99.99%
- Flow resistance: 1.5 cm H<sub>2</sub>O @ 60 L•min<sup>-1</sup>
- Low effective dead space: 45 mL
- Electrostatic media in protective membranes
- ISO 5356 Connectors: 15F/22M, 15M/22F
- Bidirectional airflow
- Disposable, for single-patient use only
- Free of latex, PVC, DEHP, BPA
- Transparent housing for visual inspection



A-M Systems strictly adheres to and manufactures products under the standards of:

- ISO 13485:2016
- FDA Quality System Regulations
- Canada Medical Device Regulations



**A-M SYSTEMS™**



SALES@A-MSYSTEMS.COM



A-MSYSTEMS.COM



TEL: +1.360.683.8300  
TEL: +1.800.426.1306



131 BUSINESS PARK LOOP  
SEQUIM, WA 98382, U.S.A.



@AMSYSTEMS1976