Ambu[®] Perfit[®] Extrication Collars



Key Benefits

- Ambu collars are fully disposable. You no longer have to worry about collars being returned to the source of origin and cross contamination issues are avoided
- Ambu collars come with a pre-shaped flip-chin piece to ensure they can be stored completely flat.
 Furthermore the flip-chin piece facilitates intubation and removal of vomit from a patient's mouth
- One-piece extrication collar simplifies application and at the same time there are no extra pieces that can get lost
- Nasal cannula holders are located on the collar body to secure the oxygen tubes (except for infant)
- All collars are radiotranslucent and come with a ventilated posterior shell for fluid drainage and a large hole for easy access. All of these features allow for patient examinations to take place safely without removing the collar

Ambu[®] Perfit[®] Extrication Collars

The Ambu Perfit collars are designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the spine during transport and routine patient care or movement.



Specifications

| Description | |
|----------------|--|
| Shell material | Polyethylene |
| Sizes | 1 (infant), 2 (paediatric), 3 (neckless), 4 (short), 5 (regular) and 6 (tall) |
| Colour | White/light grey |
| | |

Ordering specifications

| Item no. | Description |
|-------------|---|
| 264 201 000 | Perfit Infant: size 1 / cartons of 20 pcs. |
| 264 202 000 | Perfit Paediatric: size 2 / cartons of 20 pcs. |
| 264 203 000 | Perfit Neckless: size 3 / cartons of 20 pcs. |
| 264 204 000 | Perfit Short: size 4 / cartons of 20 pcs. |
| 264 205 000 | Perfit Regular: size 5 / cartons of 20 pcs. |
| 264 206 000 | Perfit Tall: size 6 / cartons of 20 pcs. |
| 264 208 000 | Carry bag for collars, fits 6 pieces / cartons of 20 pcs. |

Ambu A/S

Baltorpbakken 13 DK - 2750 Ballerup Denmark Tel. +45 7225 2000 Fax +45 7225 2053 www.ambu.com