

BURNFREE PRODUCTS

6' x 5' (183 cm x 152 cm) BurnFree Fire/Trauma Blanket



- Provides immediate pain relief by covering exposed nerve endings
- Rapidly cools and soothes the burn
- Draws the heat out and stops the burn progression
- Reduces swelling, redness, and additional skin and tissue damage
- Medical grade ingredients
- Portable, safe, and easy to use
- Provides physical protection from heat, fire, and smoke

INTL/US: UNITIZED DIMENSIONS	WEIGHT: 12.16 pounds LENGTH: 18 inches WIDTH: 9 inches HEIGHT: 7 inches
-------------------------------------	--

CASE DIMENSIONS	QUANTITY: 1 blanket WEIGHT: 13.8 pounds LENGTH: 20 inches WIDTH: 8 inches HEIGHT: 8 inches
------------------------	--

PRODUCT CODES	UPC: 702135872605 INTL/EU: N/A INTL/US: SB-7260
----------------------	---

U.S. PATENT	5,958,420
--------------------	-----------

INGREDIENTS	Purified Water, Melaleuca Alternifolia, Emulsifiers and Preservatives
--------------------	---

BurnFree solves the immediate first aid needs for all types of burns by relieving the pain, rapidly removing the heat, and cooling the burn.

The BurnFree 6' x 5' (183 cm x 152 cm) Fire/Trauma Blanket is a good size for escaping a fire. BurnFree Fire/Trauma Blankets are designed for emergency and industrial use, and to also provide physical protection from heat, fire, and smoke.

During transport to a medical facility, the blanket will soothe the victim, reducing physiological and psychological trauma. It will cool the burn, and draw the heat out to prevent further tissue and nerve damage.

The blankets are made of a pure virgin wool, woven with a unique interlinking cell construction. The special construction of the wool, allows for the blanket to hold up to 14 times its weight with BurnFree's specially formulated gel. The water-soluble gel is easily removed and will not adhere to the skin.



EXURO^{LLC}
medical[®]

BURNFREE A Division of Exuro Medical LLC
310 East 300 South, Salt Lake City, UT 84111
Ph: 801-569-9090 | Fax: 801-569-3733
info@exuromedical.com | exuromedical.com