

BURNFREE PRODUCTS

8" x 8" (20cm x 20cm) Sterile Burn & Wound Dressing



- Provides immediate pain relief by covering exposed nerve endings
- Rapidly cools and soothes the burn
- Draws the heat out and stops the burn progression
- Reduces swelling, redness, and additional skin and tissue damage
- Medical grade ingredients
- Portable, safe, and easy to use
- Sterile dressing

BurnFree solves the immediate first aid needs for all types of burns by relieving the pain, rapidly removing the heat, and cooling the burn.

BurnFree Sterile Dressings are designed for first aid use on all types of burn injuries. The dressings are constructed of a special, soft, open-cell foam, and are elastic enough to accommodate possible swelling or the need to flex joints. BurnFree dressings are non-adherent and the inert foam will not shed threads, strands, or fibers that could become embedded in the wound. The pliable foam holds 50 times its weight in gel. This means that more BurnFree gel will stay in contact with the burn.

The 8" x 8" (20cm x 20cm) BurnFree Sterile Dressing can be used for larger burns such as on the arms or face. For burns on the face, the foam can easily be cut or torn to accommodate for eyes, nose, or mouth. It is also ideal for use on larger burn injuries, such as scalding burns.

| | |
|---------------------------------|--|
| INTL/US: CASE DIMENSIONS | QUANTITY: 30 dressings WEIGHT: 10 pounds LENGTH: 13 inches WIDTH: 10 inches HEIGHT: 5 inches |
| PRODUCT CODES | UPC: 702135408088 INTL/EU: I-0808-30 INTL/US: 0808-30 |
| U.S. PATENT | 5,958,420 |
| INGREDIENTS | Purified Water, Melaleuca Alternifolia, Emulsifiers and Preservatives |



EXURO LLC
medical®

BURNFREE A Division of Exuro Medical LLC
310 East 300 South, Salt Lake City, UT 84111
Ph: 801-569-9090 | Fax: 801-569-3733
info@exuromedical.com | exuromedical.com