


## Cook Emergency Pneumothorax Set

Used for emergency relief and temporary management of suspected tension pneumothorax.



### Features and benefits

- Seldinger placement facilitates controlled, minimally invasive catheter introduction.
- The reinforced catheter is designed to prevent collapse after placement.
- The set comes with the Cook Chest Drain Valve, which is designed to permit continued re-expansion of the lung, preventing inadvertent backflow of air into the pleural space.

### Components

- One-part percutaneous entry needle
- 12 mL Monoject syringe
- Cook Chest Drain Valve
- Catheter retention disc
- Ty-Rap Cable Tie
- Vinyl connecting tube

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Needle gage	Needle Tip Configuration
G13722	C-TPTS-8.5-6.0-FSNS	8.5	6	15	Lancet
G55524	C-TPTS-8.5-8.7-FSNS	8.5	8.7	15	Lancet
<b>Pre-hospital set with ChloroPrep</b>					
G29467	C-TPTS-8.5-6.0-PHS	8.5	6	15	Lancet

Some products or part numbers may not be available in all markets. Contact your local Cook Medical sales representative or Customer Service for details.