

# Compact Triage Ribbon Kit

## MCI Optimized

Mass casualty triage has evolved and with it the tools needed to rapidly complete the task. The Compact Ribbon Kit is designed to meet the needs of MCI first responders.

### Compact Form-factor



### MOLLE / PALS Attachments



Keeping America Prepared™



### Rapid Access to Triage Supplies



#### Speed-Feed Ribbon System

Four vertically mounted 300ft ribbons feed easily but the speed-feed system prevents 'rooster-tails'.

#### Ample BLS Equipment Storage

Large user-configurable medical compartment holds critical quick-access BLS medical supplies.

#### Dual-Function Count Card

Laminated count card and included grease pens allow accurate tally in any environment. Back of card contains a quick-reference to the START process.



### Ample External Storage Points



#### Three Mesh Pockets

Front, left and right mesh pockets holds any configuration of immediate access items.

#### Bottom Bungee Sleeve

Ideal for storing additional supplies such as extra tourniquets, dressings or airway barriers.



Product No. DMS-05895