

## NYLON STRAPS

### NYLON STRAPS



Strong **NYLON** restraint straps are made to secure patients during transport.

For a custom fit, nylon straps are easily adjustable.

2" **NYLON** webbing (4,000 pound break strength)

Meets FDA requirements.



#### INSTRUCTIONS:

- 1.) Inspect strap before use.
- 2.) Position patient.
- 3.) Strap should have proper tension with no twists.
- 4.) Patient should be flat on backboard/stretchers.
- 5.) Wrap strap around patient and backboard/stretchers.
- 6.) Buckle and tighten strap to secure.

*Intended as a guide only.*

Consult your local medical agency for all protocols.

<b>1B</b> Metal Push Button Buckle 5000 Pound Break Strength	<b>5B</b> Metal D-Rings 800 Pound Break Strength	<b>9B</b> Plastic Cam Buckle 600 Pound Break Strength	<b>4E</b> Big Brass Swivel Speed Clip Ends 325 Pound Break Strength
<b>2B</b> Metal Cam Buckle 1600 Pound Break Strength	<b>6B</b> Plastic Ladderloc Buckle 550 Pound Break Strength	<b>1E</b> Loop Lock Ends	<b>5E</b> Plastic Swivel Speed Clip Ends 140 Pound Break Strength
<b>3B</b> Metal Roller Buckle 1000 Pound Break Strength	<b>7B</b> Plastic Side Release Buckle 400 Pound Break Strength	<b>2E</b> Metal Swivel Speed Clip Ends 1600 Pound Break Strength	<b>6E</b> Metal Loop Ends 495 Pound Break Strength
<b>4B</b> Metal Drop Jaw Buckle 1650 Pound Break Strength	<b>8B</b> Plastic Double Adjust Buckle 400 Pound Break Strength	<b>3E</b> Metal Fixed Swivel Speed Clip End 3800 Pound Break Strength	

<b>Black</b>	<b>Orange</b>	<b>Maroon</b>	<b>Red</b>
<b>Blue</b>	<b>Green</b>	<b>Yellow</b>	<b>Olive Drab</b>



*For use by qualified and trained personnel only.*