

NEW

Lifeline **ARM** Tactical Case

*The bag you would rather grab —
the bag you would rather stow*



PN: RAC-U1105GG

Defibtech's newest accessory* for the Lifeline ARM automated chest compression (ACC) device is a Tactical Case with a particularly unique design. The Frame and the Backboard of the Lifeline ARM encircle the custom bag. The result is a carrying case that is lightweight and takes up very little room. Visibility and access simplify basic tasks of routine maintenance.

Check Battery Pack Status Indicator – Access Input Jack to Recharge

- Clear window on top of the bag provides visibility of the Battery Pack status indicator on the User Control Panel
- To Recharge the device while it is in the bag, access the input jack of the Compression Module via a special flap opening on the front side of the bag

Bag is made of padded, durable materials, including reinforced seams throughout. Inside the bag, an interior pocket secures and protects the Compression Module, with additional pockets and space for storing accessories and spares. Comfortable carrying options include a center top handle, rear handle, and padded backpack straps. Overall, the durability, accessibility, and comfort of this carrying case are just a few of the things that set this bag apart.

Quickly know you are ready to go



Lightweight, compact, and you can check Battery Pack charge level, or recharge the Battery Pack, through the bag

It's in the bag:

- Red exterior with hi-visibility retro-reflective stripes for low light conditions
- Brightly colored interior makes it easier to see stored accessories and spares
- Suitable for simple decontamination
- Weight 1.25 lbs. (0.57 kg); 22.5" L x 12" H x 9.5" D (case only)

**Optional accessory with 1-year warranty*

Life-Assist
Helping Heroes Save Lives

defibtech
Life-Saving Technology Within Reach

Defibtech, LLC | 741 Boston Post Road, Guilford, CT 06437 USA
www.defibtech.com | 1-866-DEFIB-4U (1-866-333-4248) or 1-203-453-4507

RAC-U1705EN-BA
ELECTRONIC DISTRIBUTION