










## ECG Electrode Standardization Chart

Selecting the correct electrode can be challenging. CONMED has designed its electrode portfolio to be all-inclusive and meet the most common patient and clinician needs, while also offering an array of specialty electrode products to meet the needs of most hospital department areas.

**One Supplier for All Your Electrode Needs**

Patient Type	Adult				Tabs	Specialty	Pediatric/Neonatal		
									
<b>Ordering Number Name</b>	1870   1870C Positrac	2700   2710 Cleartrac 2	1700   1700C Cleartrac	3700   3700C Totaltrac	1915 Fastrac	2330   2330C Softtrac	1720   1720C Cleartrode	1741   1741C Neotrode	2331   2331C Softtrac Pre-attached Leads
<b>Adhesion Strength</b>	Diaphoretic	Diaphoretic	General Purpose	General Purpose	Repositionable	Gentle	Gentle	Gentle	Gentle
<b>Basepad</b>	Foam	Tape	Tape	Foam	Mylar	Cloth	Tape	Tape	Cloth
<b>Gel Type</b>	Solid Gel	Solid Gel	Solid Gel	Solid Gel	Solid Gel	Repositionable	Repositionable	Repositionable	Repositionable
<b>MR Safe</b>	1870 - No 1870C - Yes	2700 - Yes 2710 - No	1700 - No 1700C - Yes	3700 - No 3700C - Yes	No	2330 - No 2330C - Yes	1720 - No 1720C - Yes	No	No
<b>Radiotranslucent</b>	1870 - No 1870C - Yes	2700 - Yes 2710 - No	1700 - No 1700C - Yes	3700 - No 3700C - Yes	Yes	2330 - No 2330C - Yes	1720 - No 1720C - Yes	1741 - No 1741C - Yes	2331 - No 2331C - Yes
<b>Pouch Quantity</b>	1870- 1, 3, 5, 30, 50 1870C- 3, 4, 5, 50	1, 3, 5, 30	1700- 1, 3, 4, 5, 30, 50 1700C- 3, 5, 10, 50	1, 3, 5, 50	10, 100	3	3	3	3