

Sta-Blok™ Head Immobilizer



Laerdal Solutions for EMS

www.laerdal.com



Laerdal[®]
helping save lives

Sta-Blok™ Head Immobilizer


Foam immobilizers are comfortable for the patient and easy for the provider to use. As part of Laerdal's total spinal solution, we are pleased to offer the Sta-Blok Head Immobilizer. Compact, adjustable, effective, and durable yet disposable, Sta-Blok is carefully constructed of soft, closed-cell foam and can fit virtually any patient on any spineboard. (Of course, we prefer the BaXstrap™ Spineboard, one of the lightest, more comfortable boards available.)

- **Single-use disposable devices** means no equipment recovery costs, require no clean-up, minimize contamination and can be used in mass casualty situations
- **Quick and easy application** minimizes time spent on the scene

- **Easy access to ears and head** without removing the device for ongoing patient assessment
- **Radiolucent:** allows patient to remain in device for X-Ray evaluations, minimizing patient movement
- **“Severe Weather/Helmet” strap** provides extra stability and can be used in any type of weather condition
- **Adjustable, standard Velcro® padded strap** is easy to re-adjust, allows for quick access and eliminates tape on the patient's eyebrows
- **Latex Free**

Sta-Blok Features:

- Velcro padded strap – no tape on eyebrows
- Constructed of soft closed-cell foam
- Easy to adjust



Universal Base Features:

- Helps elevate spine to neutral position
- Provides patient comfort

700-00001 Sta-Blok Head Immobilizer (pkg. 30)

“Issues in Spinal Care” Training Curriculum CD 2005 Edition

The classic spinal care training curriculum is now revised and updated!

The 2005 Edition presents concepts and techniques of spinal motion restriction in four independent modules: adult, pediatric, aquatic, and spinal field rule-out. Most of the updates are in the aquatics section. Independently and expertly reviewed by the Minnesota State Colleges and Universities Fire/EMS/Safety Center and the physicians listed on their Medical Advisory Committee.



For more information, visit us at www.laerdal.com. To order call 877-LAERDAL (523-7325) or contact your local Laerdal Representative.



Laerdal®
helping save lives