

EMESIS BAG

#3933

*Convenient personal solution that
can be easily placed anywhere
in the facility*

FEATURES:

Convenience: Wide opening is fully supported by a rigid plastic ring to help ensure unobstructed access to the containment bag while providing a circular grasping surface for secure grip and proper positioning around the mouth.

Capacity: Deep cylindrical shape holds in excess of 1000 cc/ml (40 fl. oz.) for managing vomit to help prevent spills, and costly clean ups. Durable polyethylene construction is lightweight and tear-resistant.

Minimize Exposure: Rigid ring support channels fluid into the containment bag to help reduce spillage and manage exposure.

Secure Disposal: Twist, close, and hook bag into plastic ring support notch for easy, one-step disposal of contents and emesis bag. Odors are contained during transport.

Easy Dispensing: Wall Hanging Dispenser (#3933D) holds folded containment bags securely and facilitates individual dispensing when needed. Dispensers provided with minimum annual purchase commitment of Emesis Bags, or can be purchased separately. Contact your Medegen Medical Products sales representative for details.

Areas of Use: Install virtually anywhere patient services are provided.



MEDEGEN
Medical Products
An Inteplast Group Company

209 Medegen Drive, Gallaway, TN 38036
Toll Free 800-511-6298 • Customer Service Fax 631-714-6311
E-mail cservice@medegenmed.com



1 Always within easy reach when used with the Wall Hanging Dispenser (#3933D).




2 Rigid plastic ring provides a wide opening and secure grasping surface.



3 Twist used Emesis bag to safely contain fluids.



4 Dispose after hooking twisted bag into the support ring notch.



Emesis Bag (#3933) Specifications

- Dimensions: 6.5" X 15" overall, containment bag length 14"
- Film: LDPE, transparent blue tint allows content visualization
- Graduations: To 1000cc/ml and 40 fl. oz.
- Case Quantity: 100 each/case, 25 bags per sleeve (4)
- Case Weight: 3.8 lbs.
- Latex Free

Wall Hanging Dispenser (#3933D)

- Ideal For:
- Emergency Room
- Oncology
- Radiology
- PACU Recovery
- Labor and Delivery
- Pediatrics
- General Medical
- Waiting Areas

Dispenser for use with Medegen Medical Products part# 3933 Emesis Bag refills only



For more information, please contact your Medegen Medical Products sales representative at 1-800-511-6298. To view our products online, please visit us at www.medegenmed.com or directly with your device by scanning the QR code below.

Medegen Medical Products develops, manufactures and distributes a comprehensive portfolio of high quality disposable medical and surgical products to the healthcare community. Medegen is a leading provider of disposable products, patient bedside products, containment systems for medical waste and protective apparel.



Medegen Medical Products
209 Medegen Drive
Galloway, TN 38036

CUSTOMER SERVICE 800-511-6298
CUSTOMER SERVICE FAX 631-714-6311
CUST. SERVICE EMAIL CService@medegenmed.com

INTL CUST. SERVICE IntCService@medegenmed.com
WEB www.medegenmed.com

© 2015 Medegen Medical Products
LIT# 0037

Life-Assist
Helping Heroes Save Lives

INTEPLAST GROUP
Inteplast Healthcare