

MICRO-TOUCH[®] NITRAFREE[™]

H i g h T o u c h

Provides comfortable and reliable protection while supporting breast cancer research

The MICRO-TOUCH[®] NITRAFREE[™] Nitrile examination gloves are made from a durable and soft nitrile formulation manufactured without chemical accelerators, helping to protect users from Type IV chemical allergies. Provides comfortable and reliable protection while supporting breast cancer research.

Recommended For

- Tested for use with chemotherapy drugs
- Standard patient examinations
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] NITRAFREE[™]

High Touch

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Pink	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.14 / 5.5	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Silicone Free	Not Tested	Sterile	No
Antistatic	Not Tested	Double Gloving Recommendation	Yes, as underglove
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034510	6034511	6034512	6034513	6034514

PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, EN 420:2003 + A1:2009, EN 421:2010, EN 455 1-4, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600, ISO 9001, ISO 13485, ISO 10993-10:2002

Vulcanization Chemical Accelerators

None	
------	--

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email: insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking.
© 2021 Ansell Limited. All Rights Reserved.

ansell.com