

O2-MAX™

Disposable CPAP Device


Continuous positive airway pressure is becoming the standard of care for patients suffering from respiratory distress. The O2-MAX system can deliver up to 140 LPM of flow to the patient, with 30% FiO₂. Customize your FiO₂ to increase from 30, 60, or 90% with the Trio adapter. Adjustable PEEP options include a 3-SET™ (5/7.5/10cm H₂O) or 5-SET (5/7.5/10/12.5/15cm H₂O) to suit a variety of patient needs. Additionally, this completely disposable system offers integrated nebulization while maintaining a consistent flow to the patient.



28 Minutes of Run Time **

Adjustable PEEP (3-SET™ & 5-SET™)

Filter Efficiency:
99.9998% Bacterial Efficiency
99.9998% Viral Efficiency



**D-Cylinder Tank with 30% FiO₂

O2-MAX™ Fixed System

Part No.	PEEP Setting	Size	FiO ₂	Circuit	Nebulizer	Quik-Connect	Case/Qty.
313-7556XN-1	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30%	Expandable	Yes	Ohmeda	10
313-7555XN-1	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30%	Expandable	Yes	Ohmeda	10
313-7557XN-1	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30%	Expandable	Yes	Ohmeda	10
313-7554XN-1	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30%	Expandable	Yes	Ohmeda	10
313-7556XN	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30%	Expandable	Yes	-	10
313-7555XN	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30%	Expandable	Yes	-	10
313-7557XN	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30%	Expandable	Yes	-	10
313-7554XN	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30%	Expandable	Yes	-	10
313-7556X	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30%	Expandable	-	-	10
313-7555X	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30%	Expandable	-	-	10
313-7557X	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30%	Expandable	-	-	10
313-7554X	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30%	Expandable	-	-	10

O2-MAX™ Trio System

Part No.	PEEP Setting	Size	FiO ₂	Circuit	Nebulizer	Quik-Connect	Case/Qty.
313-8007XN-1	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Expandable	Yes	Ohmeda	10
313-8007XN-6	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Expandable	Yes	DIN (German)	10
313-8017XN-1	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Expandable	Yes	Ohmeda	10
313-8017XN-6	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Expandable	Yes	DIN (German)	10
313-8007N-5	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	Yes	American Gas	10
313-8008N-5	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	Yes	American Gas	10
313-8008N-6	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	Yes	DIN (German)	10
313-8017N-5	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	Yes	American Gas	10
313-8018N-5	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	Yes	American Gas	10
313-8018N-6	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	Yes	DIN (German)	10
313-8028N-13	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	Yes	AFNOR & DIN (German)	10
313-8038N-13	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	Yes	AFNOR & DIN (German)	10
313-8007-5	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	-	American Gas	10
313-8017-5	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	-	American Gas	10
313-8007XN	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Expandable	Yes	-	10
313-8008XN	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30/60/90%	Expandable	Yes	-	10
313-8017XN	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Expandable	Yes	-	10
313-8018XN	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30/60/90%	Expandable	Yes	-	10
313-8007X	3-SET™ (5/7.5/10cm H ₂ O)	Adult Medium	30/60/90%	Expandable	-	-	10
313-8017X	3-SET™ (5/7.5/10cm H ₂ O)	Adult Large	30/60/90%	Expandable	-	-	10
313-8008	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Medium	30/60/90%	Smooth Bore	-	-	10
313-8018	5-SET™ (5/7.5/10/12.5/15cm H ₂ O)	Adult Large	30/60/90%	Smooth Bore	-	-	10

O2-MAX™ Accessories

Part No.	Description	Size	FiO ₂	Case/Qty.
313-7028	BiTrac ED Mask with Standard Headgear	Adult Small	-	10
313-7029	BiTrac ED Mask with Standard Headgear	Adult Medium	-	10
313-7030	BiTrac ED Mask with Standard Headgear	Adult Large	-	10
313-7031	BiTrac ED Mask with Standard Headgear	Adult Extra Large	-	10
313-8000	O2-MAX Trio Adapter	N/A	30/60/90%	10