



TACMED
SOLUTIONS
ELEVATE YOUR RESPONSE



 **TACMED**
MEDICAL

OLAES[®] MODULAR BANDAGE

THE MOST VERSATILE TRAUMA BANDAGE ON THE MARKET.

PRODUCT SPECIFICATIONS

PRODUCT #	4" Version - OAL-4 4" Flat Packed - OALFP-4 6" Version - OAL-6
DIMENSIONS	4" x 4" Flat Packed - 4"W x 7"L x 2.5"D 10.2cm W x 17.8cm L x 6.4cm D 6" Version - 4"W x 7.75"L x 2.5"D 10.2cm W x 19.7cm L x 6.4cm D
WEIGHT	4" x 4" Flat Packed - 3oz / 85g 6" Version - 3.7oz / 105g
NSN	4" - 6510-01-558-3342 4" Flat Packed - 6510-01-631-7893 6" - 6510-01-558-3443



INTENDED USE

For treatment of moderate hemorrhage.

Life-Assist
Helping Heroes Save Lives 

OLAES[®] MODULAR BANDAGE



The **OLAES[®] Modular Bandage** is advantageous for all levels of care providers and is the most comprehensive feature-packed bandage available. The OLAES[®] provides the ability to cover minor wounds, apply truly effective pressure dressings, apply occlusive dressings, and pack severe wounds. All of the advanced capabilities of the bandage are tucked neatly away and can be utilized if required or simply ignored.

PRODUCT FEATURES

- 1 HOOK & LOOP CONTROL STRIPS**
Built into the elastic wrap to prevent accidental unrolling during application.
- 2 TRUE PRESSURE CUP**
Generates focused pressure on the wound site to assist occluding damaged vessels. It can also serve as a rigid eye shield for ocular injuries per CoTCCC guidelines.
- 3 OCCLUSIVE SHEET**
Plastic, occlusive sheet located behind the wound pad for increased treatment versatility.
- 4 STERILE GAUZE**
Stocked with 3 meters of sterile, zig-zag packed gauze inside the wound pad, ideal for packing wounds.

- 5 WOUND PAD**
Covers the injury site to assist with sustained bleeding control and keep the wound clear of debris. Its non-adherent compartment keeps gauze tucked neatly inside for rapid access during application.
- 6 SECURING CLIP**
Locks the bandage into place to free up the user's hands to treat other injuries or casualties.

BANDAGE OVERALL LENGTH

4" - 48.0"±2"/ 121.92cm±5.08cm
4" Flat Packed - 48.0"±2"/ 121.92cm±5.08cm
6" - 48.0"±2"/ 121.92cm±5.08cm

ADDITIONAL FEATURES

SIMPLER SOLUTION

During high stress situations, an individual's fine motor skills can be affected, so the OLAES[®] only involves gross motor skills to offer a significantly simpler and more effective bandage solution for applications.

SHELF LIFE

6 Years from Date of Sterilization

MATERIALS LIST

ABS, Cotton (8-Ply), Elastic (Latex Free), Nylon, Polyester (Non-Woven), Polyethylene (PE), Polyethylene Terephthalate (PET), Rayon (Thermo-Bonded Filler).