

Malleability

Oxygenation*

Soft Tip



S-Guide

Malleable Intubating Guide
for Difficult Airway Management

S-Guide – 3 benefits in 1



Application video

The S-Guide is an intubation aid for difficult airway management and facilitates orotracheal intubation, especially in patients where visualisation of the glottis is inadequate. Due to its **malleability**, it is suitable for direct and indirect laryngoscopy, in particular for video laryngoscopes with non channelled blades. With the O₂ connectors, **oxygenation*** is possible throughout the whole intubating process, not impairing oxygenation even when the S-Guide is bent. For patient safety the S-Guide is provided with a coloured, preformed and atraumatic **soft tip**.

Malleability

- Specific design of material combines:
 - Malleability, to preform the tracheal tube
 - Rigidity, to manoeuvre the tracheal tube
- Special black marking on the sizes 11 Fr and 15 Fr shows ideal position to preform the tracheal tube (hockey-stick marking)



Oxygenation*

- For sizes 11 Fr and 15 Fr
- Three openings at the distal end of the S-Guide to ensure appropriate oxygen flow
- 2 different connectors:
 - O₂ - 15 mm Connector and O₂ Connector



Soft Tip

- Preformed, soft and atraumatic tip
- Coloured tip:
 - For better visualisation within the glottic structures
 - For size coding



*Oxygenation is possible with the sizes 11 Fr and 15 Fr. Oxygenation is not possible with the size 8 Fr.

Order Information

S-Guide

| Size | Length | Tracheal tube I.D. | REF | Box |
|-------|--------|--------------------|--------------------|-----|
| 8 Fr | 40 cm | ≥ 3.0 mm | 33-92-400-1 | 5 |
| 11 Fr | 65 cm | ≥ 4.5 mm | 33-91-650-1 | 5 |
| 15 Fr | 65 cm | ≥ 6.0 mm | 33-90-650-1 | 5 |

Product Specifications

- Shelf life of 5 years
- Sterile packed
- Single use product

This product is manufactured without the use of natural rubber latex, unless otherwise specified.
This product does not contain phthalates which require labelling according to CLP Regulation (EC) 1272/2008.

VBM Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany

Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.de

Connect with us:



CE 0123

PRO019_2.0_EN