



# ChitoSAM™

[kahy-tuh-sam]

# COVER. CLOT. CONTROL.

## STOP BLEEDING RAPIDLY


ChitoSAM™ 100 is a high-performance hemostatic, non-woven chitosan dressing spun directly from chitosan derived crustaceans or snow crab shells. Designed to stop lethal bleeding rapidly, its ease of use is extremely effective. Works independently from the body's normal clotting processes.



## FEATURED TECHNOLOGY

### TEAR ANYWHERE PACKAGING

Innovative package design allows for quick access to ChitoSAM™ 100, tearable from any side.



**ENGINEERED TO PRESERVE LIFE**

- 1 PERFORMANCE GAUZE**  
ChitoSAM<sup>™</sup> 100 has a hemostatic performance equivalent to products impregnated or coated with chitosan hemostatic agents.
- 2 CHITOSAN HEMOSTATIC AGENT**  
Unlike most hemostatic dressings that are only coated or impregnated with chitosan, ChitoSAM<sup>™</sup> 100 is made from non-woven chitosan gauze, spun directly from chitosan derived from crustaceans or snow crab shells.
- 3 TEAR ANYWHERE PACKAGING**  
Innovative package design allows for quick access to ChitoSAM<sup>™</sup> 100, tearable from any side.
- 4 Z-FOLD**  
Available for large wounds in 3in x 4ft and 3in x 6ft.
- 5 SINGLE PLY**  
Smaller 4in x 4in single ply ChitoSAM is able to be placed on wounds in addition to trauma gauzes.



\*CE versions shown

**PRODUCT SPECIFICATIONS**

**CHITOSAM<sup>™</sup> 100 MODELS AVAILABLE (\*CE)**

ChitoSAM 100 4in x 4in (10cm x 10cm) (CTS101-44)  
ChitoSAM 100 3in x 6ft (7.6cm x 1.83m) (CTS101-06)

**CHITOSAM<sup>™</sup> 100 MODELS AVAILABLE (Non-CE)**

ChitoSAM 100 4in x 4in (10cm x 10cm) (CTS100-04)  
ChitoSAM 100 3in x 4ft (7.6cm x 1.22m) (CTS100-44)  
ChitoSAM 100 3in x 6ft (7.6cm x 1.83m) (CTS100-06)

**WARNING**

Do not place over eyes.

**GLOBAL SALES**

U.S. Toll Free: 800.818.4726  
U.S. Local: 503.639.5474  
U.S. Fax: 503.639.5425  
sammedical.com

<b>EC REP</b>	<b>CE</b> 0050	<b>2</b>	<b>R ONLY</b>	<b>STERILE</b>	<b>R</b>
Emergo Europe, Prinsessegracht 20, 2514 AP, The Hague, The Netherlands				SAM Medical Products <sup>®</sup> 27350 SW 95th Ave, Ste 3038 Wilsonville, OR 97070 USA	

SAM<sup>®</sup> is a registered trademark in the United States and other countries.  
©2017 SAM Medical Products. All rights reserved.